English

## Disposable Extended Monopolar Stimulating Probe

### Intended Applications

IONM [Intraoperative Neurophysiological Monitoring].

### Caution

Monopolar Stimulating Probes are for professional use only and should only be used in compliance with accepted Industry standards. Rhythmink International, LLC is not responsible for injury, infection or other damage resulting from the use or misuse of this product. Product is sterile only if packaging is unopened. Do not use near high energy radio frequencies.

### Intended Use

Rhythmink International Monopolar Stimulating Instrument is indicated for stimulation of cranial and peripheral motor nerves for location and identification during surgery, including spinal nerve roots. The Rhythmink International Monopolar Stimulating Instrument is sterile and for single use only.

### Instructions for Use

Remove protective cover from Probe. Insert the connector end into the stimulation device socket. This product is single-patient use only. Discard Probe after use.

Rhythmink International, LLC  
1140 First Street South  
Columbia, SC, USA 29209-3540  
+1.866.633.3754 [toll-free]  
+1.803.252.1222  
+1.803.252.1111 [fax]  
sales@rhythmink.com  
Rhythmink.com

**DO NOT BEND**

**STERILE EO**

**EC REP**  
European Authorized  
Representative  
Mdi Europa GmbH  
Langenhagener Str.71  
D-30855 Hannover-Langenhagen

- Reuse of this device may cause cross contamination and has NOT been validated for multiple uses.
- Re-sterilization of this device may cause unknown malfunctions. EO re-sterilization has not been validated for this device and may contribute to excessive EO gas residuals.

To view a list of symbol definitions found on packaging and instructions for use, please visit [Rhythmink.com/symbols](http://Rhythmink.com/symbols).

Rhythmink® is a registered trademark of Rhythmink International, LLC.

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Rhythmink® est une marque déposée de Rhythmink International LLC
Rhythmink® ist eine eingetragene Marke von Rhythmink International, LLC.