English

Disposable PressOn™ Electrode Headset

Intended Use

The PressOn™ Electrode is intended for use in the recording of the Electroencephalogram (EEG), the evoked potential (EP), or as a ground and reference in an EEG or EP recording. This device is provided sterile for Single Patient Use Only.

Caution

Federal [USA] law restricts this device to sale by or on the order of a physician and it should only be used in compliance with accepted industry standards. RhythmLink International, LLC is not responsible for injury, infection or other damage resulting from the use or misuse of this product. Product is sterile unless packaging is opened or damaged.

Packaging

A sterile pouch contains electrodes with positioning headset. Extension cables are packaged separately in individual sterile pouches.

General Information

- PressOn electrodes are supplied sterile, intended for single patient use and are disposable.
- The PressOn EEG leads can be connected to any commercially available EEG equipment.
- PressOn electrodes should only be used by trained or skilled personnel under the direction of a physician.
- PressOn electrodes are approved for prolonged usage.
- Follow hospital guidelines for application and dispensing of the PressOn electrode.

Warnings

- Electrodes are only approved for recording.
- Electrodes should be applied only to normal, healthy skin.
- The skin, shape and type of electrodes may affect the safety and effectiveness of electrical recording.
- The electrical performance characteristics of electrodes may affect the safety and effectiveness of electrical recording.
- You should contact the manufacturer of the electrical recording device if you do not know if the electrodes can be used with the recording device.
- The entire device is MR Unsafe.

RhythmLink International, LLC
1140 First Street South
Columbia, SC, USA 29209-3540
+1.866.633.3754 [toll-free]
+1.803.252.1222
+1.803.252.1111 [fax]
sales@rhythmLink.com
RhythmLink.com

FOR SINGLE USE ONLY

STERILE

EO

Pre-Sterilized Disposable

Rx Only

To view a list of symbol definitions found on packaging and instructions for use, please visit [RhythmLink.com/symbols](https://rhythmLink.com/symbols).

Reuse of this device may cause cross contamination and has NOT been validated for multiple uses.


Re-sterilization of this device may cause unknown malfunctions. EO re-sterilization has not been validated for this device and may contribute to excessive EO gas residuals.

RhythmLink® is a registered trademark of RhythmLink International, LLC. PressOn™ is a trademark of Persyst Development Corporation.

Precautions

The long term effects of cutaneous electrodes for electrical recordings are unknown. Keep electrodes out of the reach of children. Use caution if electrodes are applied over areas of skin that lack normal sensation.

Instructions for Use

1. Open the sterile pouch and remove the headset from the pouch.
- 2. Position the headset on the patient using the chin strap to secure. Specific applicators indicate the front of the headset and should rest on the patient's forehead.
3. Adjust the headset as needed so it rests evenly and securely on the patient. The applicators should be oriented so that the electrode end is resting on the skin and the applicator's button is facing up. Take care not to deploy electrodes until ready.
4. If necessary, part patient's hair at the electrode placement site to ensure the PressOn is touching the skin prior to launch.
- 5. Stabilize each applicator with your thumb and middle finger and press the plunger downward with your pointer finger toward the skin with light to moderate pressure to apply the electrode. If properly applied the teeth should be equally embedded under the skin.
6. When all electrodes are deployed, release the chin strap and gently remove the positioning headset using caution not to dislodge any electrode.
7. Dispose of positioning headset according to standard biomaterial handling protocols.
- 8. Remove extension cables from the pouch and match the numbers on the extenders to the patient array connector to ensure accuracy. Next, plug the electrodes into any commercial EEG equipment per manufacturer's instructions.
9. During prolonged monitoring, check the electrode insertion site regularly for any signs of infection, reddening or discharge.
- 10. To remove a PressOn electrode, grasp the electrode's wire within 1-2 centimeters of the electrode site, then pull quickly upward away from the skin. Pull in a direction perpendicular to the skin until the electrode pulls free. Discard all used electrodes per hospital guidelines.
11. Minor capillary bleeding may be present once an electrode has been removed. Should this occur, press a sterile cotton pad on the site and hold until bleeding has stopped or follow hospital guidelines for care.
12. Dispose of electrodes in a medical sharps container according to standard biomaterial handling protocols.

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