



RhythmLink International, LLC
1140 First Street South
Columbia, SC 29209
RhythmLink.com

866.633.3754 (toll-free)
803.252.1222
803.252.1111 (fax)
sales@rhythmLink.com

Disposable Single Sticky Pad™ Surface Electrodes

Materials

- Pad - Conductive Vinyl layer coated with Cotton Micropore
- Pad Backing - Mylar
- Stimulating Hydrogel - AG702 Stimulating
- Recording Hydrogel - AG603 Sensing
- Ground Hydrogel - AG603 Sensing
- Leadwire - Tinned Copper with PVC Jacket
- Connector - DIN 42 802
- Sterilized - No

Product Specifications

- Stimulating and Recording Pad Sizes - 15 x 20 and 20 x 25 mm
- Ground Pad Size - 45 x 35 mm
- Leadwire Lengths - 1.0, 1.5 or 2.5 m

Environment

- Latex-Free Electrode
- PVC-Free Packaging

Packaging

- Pouch - Metalized Foil
- Pouch Size - 130 x 150 mm
- Boxes - Cardboard

Additional Benefits

- Stimulating pads feature a specially formulated stimulating hydrogel that minimizes stimulation risks
- Recording pads feature a specially formulated recording hydrogel that optimizes signal recording
- Identifying letters distinguish between pads
- All Electrodes feature molded safety connectors with ample built-in strain relief

RhythmLink's Disposable Sticky Pad Surface Electrodes are cleared for both stimulating and recording and available in three varieties: Stimulating, Recording and Ground. These versatile, robust pads feature an aggressive hydrogel which keeps the Electrodes in place and enhances physiologic recordings or stimulations.

Visit RhythmLink.com/stickypad for additional information.

Directions for Use - Select appropriate Sticky Pad Electrode type and size. Remove excess hair, oils and dirt. Prep application site, remove Mylar backing from Electrode and place on application site. Apply pressure to center of Electrode and move to edges. Remove by pulling directly on Electrode. **DO NOT PULL ON THE CABLE!** Remove remaining hydrogel with clean, soapy water.

Cleaning & Disinfecting Instructions - This product is single-patient use only. Discard Electrode after use.

Other Sticky Pad Varieties include - Paired Stimulating and Paired Recording

