



**RhythmLink International, LLC** 866.633.3754 (toll-free)  
1140 First Street South 803.252.1222  
Columbia, SC 29209 803.227.1015 (fax)  
[RhythmLink.com](http://RhythmLink.com) [sales@rhythmLink.com](mailto:sales@rhythmLink.com)

## Disposable Alligator Clips - Small

### Materials

- Clip - Stainless Steel
- Hub - Polyolefin
- Hub Diameter - 1.25 mm
- Leadwire - Tinned Copper with PVC Jacket
- Connector - DIN 42 802
- Sterilized - Yes, Ethylene Oxide

### Product Specifications

- Clip Size - 43 mm
- Leadwire Length - 2.5 m

### Environment

- Latex-Free Electrode
- PVC-Free Packaging

### Packaging

- Pouch (*Transparent Layer*) - Nylon/Polyethylene
- Pouch - Tyvek
- Pouch Size - 86 x 160 mm
- Box - Paperboard

RhythmLink's Alligator Clips are a great accessory for general surgical instrumentation. A small clip, this product is 100% disposable and intended for use in the recording of intraoperative neuromonitoring.

Visit [RhythmLink.com/alligatorclips](http://RhythmLink.com/alligatorclips) for additional information.

**Directions for Use** - The Alligator Clip is intended for use in the recording of IONM as an accessory for general surgical instrumentation.

**Cleaning & Disinfecting Instructions** - This product is single-patient use only. Discard clip after use.

