



Rhythmlink International, LLC 866.633.3754 (toll-free)
1140 First Street South 803.252.1222
Columbia, SC 29209 803.252.1111 (fax)
Rhythmlink.com sales@rhythmlink.com

Disposable Webbed™ EEG Electrodes*

Materials

- Cup - ABS Plastic Coated with Ag/AgCl
- Leadwire - Tinned Copper with PVC Jacket
- Connector - DIN 42 802
- Sterilized - No

Product Specifications

- Plating - Ag/AgCl
- Cup Size - 10 mm
- Leadwire Lengths - 1.0, 1.5, 2.0 and 2.5 m

Environment

- Latex-Free Electrode
- PVC-Free Packaging

Packaging

- Pouch - Polyethylene
- Pouch Size - 120 x 182 mm

Additional Benefits

- A molded safety connector with ample built-in strain relief
- Pre-attached leadwire
- Can be used with any gel, paste, collodion or other similar adhesive
- Eliminate cross contamination

Rhythmlink's Disposable Webbed EEG Electrodes are designed to increase patient comfort while providing quality recordings that can last through both night and day studies and for multiple days. As a disposable electrode, remove the risk of cross contamination by using it once and throw it away.

Visit Rhythmlink.com/eeg for additional information.

Directions for Use - Clean application site. Apply Electrode using conductive gel, paste, collodion or other similar adhesive. When finished, remove Electrode and clean application site.

Cleaning & Disinfecting Instructions - This product is single-patient use only. Discard Electrode after use.

* Patent No. D 644, 738

