



RhythmLink®

PRODUCT CATALOG

Neurodiagnostic Solutions
EEG, EMG, IONM and Sleep Products



OUR MISSION & VISION

To enhance patient care worldwide by transforming the medical device technology that links patients to machines.

Rhythmink International, LLC designs, manufactures and distributes medical devices and provides custom packaging, private labeling, custom products and contract manufacturing to its customers. Rhythmink is recognized as a leader within its field at providing the important physical connection between patients and the diagnostic equipment to record or elicit neurophysiologic biopotentials.

Originally founded by neurodiagnostic technicians and engineers in 2002, Rhythmink strives to provide continuous innovation and superior quality in all of its products. Based in Columbia, SC, Rhythmink's advancements and improvements in technology, business development and corporate branding have brought national and international recognition.

We deliver peace of mind by providing superior products, consistent availability, continuous innovation and the highest level of customer service in the industry without wavering in our commitment to fair pricing.

To this end, we are Rhythmink International, LLC.

REQUESTING SAMPLES

Request samples today at Rhythmink.com, sales@rhythmink.com or 866.633.3754.

STANDING ORDERS

Remove the time spent calling each month to generate the same order by setting up a Standing Order. Our Sales Team will work with you to create your recurring order which will then automatically ship to you by the 20th of each month.

COMPETITIVE PRICING

Rhythmink offers competitive pricing based on volume with our promise to deliver your products promptly. We are happy to work with your service to meet your needs.

SPECIAL ORDER & CUSTOM PRODUCTS

Special Order and Custom Products may have an additional lead time due to the nature of custom manufacturing. These products are not part of the sample program. Contact our Sales Team for more information at sales@rhythmink.com.

DON'T SEE WHAT YOU'RE LOOKING FOR?

Rhythmink specializes in transforming ideas into real life solutions. Whether it's modifying an existing product or creating a brand new one, partner with Rhythmink to create custom solutions specific to your needs. Learn more about our custom options on page 46.

TABLE OF CONTENTS

EEG Products

MR Conditional/CT Electrodes	
Quick Connect System™	5
Individual Connectors	20
Sticky Pad Surface Electrodes™	23
Artifact Free CT Electrodes	24
RLI EEG Template	25
RLI BrainHealth™ Headset	25
PressOn™ Electrodes	26
Standard EEG Electrodes	26

EMG/NCV Products

Monopolar Needles	32
Concentric Needles	33
NCV Sticky Pad Surface Electrodes	34

IONM Products

Stimulating Probes	36
Subdermal Needles	
Single Needles	38
Paired Needles	39
Twisted Paired Needles	40
5 & 6 Twisted Needles	40
Bent Needles	40
Corkscrew Needles	41
Fixed Distance Guardian Needle	41
Guardian Needle	42
Standard Sticky Pad Surface Electrodes	43
Alligator Clips	44
ET Tube Electrodes	44
SightSaver™	45
Gels & Pastes	45
Product Listing by Part Number	47
Terms and Conditions	49

RhythmLink - EEG Products



RHYTHMLINK EEG ELECTRODES

RHYTHMLINK ELECTRODE STYLES

Deep Cup™	Slim Cup™	Webb™	Small Webb™	PressOn™
- Traditional cup - Designed to hold ample conductive material - Incorporates a broader rim design to assist with the protection of skin health	- Lower profile dome to assist with patient comfort - Incorporates a broader rim design to assist with the protection of skin health	- Flat electrode design for increased patient comfort - Greater surface area designed with skin health at the forefront - Perfect for areas on the patient where skin health or skin injury is at the highest risk	- Our same high quality flat electrode design in a smaller size - At just 6 mm wide, it's perfect any time you need a smaller electrode - Perfect for delicate skin or areas on the patient where skin health or skin injury is a concern	- Designed to be placed quickly for fast, quality recordings - Designed to utilize medical staff already at bedside to collect real time data and provide the best patient care possible
				
- Standard - CT - MR Conditional	- Standard - CT - MR Conditional	- Standard - CT - MR Conditional	- Standard - CT - MR Conditional	- Standard - CT - MR Conditional

STANDARD, CT OR MR CONDITIONAL

1. Standard EEG Electrodes - Disposable electrodes designed to provide quality recordings in a variety of styles and can be used for most CT studies.

2. CT Invisa-Electrode™ - Our CT design is an artifact free disposable electrode that can be left on your patient during a CT scan.


3. MR Conditional/CT - FDA cleared MR Conditional for 0.2, 0.5, 1.0, 1.5 & 3T. Can safely remain on your patient during an MRI or CT scan.

Individual Connectors - Offer longer lengths and more flexibility for anatomical needs.


Quick Connect System™ - Offer fewer connection points for quick and efficient disconnection and reconnection of electrodes when placing, plus improved cable management.

HOW TO USE RHYTHMLINK'S QUICK CONNECT SYSTEM

When a patient is referred to MRI simply disconnect by pressing down gently on the yellow tab and pull apart.



When the patient returns from the radiology suite, use the number labeling to reconnect. Matching number to number line up the connectors and insert the tab. You will hear a click which indicates the connector is locked in place.



See our full Quick Connect User Video online at www.rhythmink.com/QC

WHY RHYTHMLINK


RhythmLink is the Market Leader in MR Conditional EEG Electrodes, having set the standard over ten years ago with our Individual Connectors. We've now manufactured over 30 million MR Conditional EEG Electrodes and provide the widest offering of already available SKUs, plus customization options to match your protocols.



**SAVE TIME AND PROVIDE
MORE EFFICIENT CARE
COMPARED TO STANDARD EEG**

WAYS TO CUSTOMIZE

If you don't see a configuration that's perfect for you, let RhythmLink make a fully custom set of electrodes to match your workflow. Fully customize each of the following components:



**Ready to design your custom configuration?
Contact sales@rhythmLink.com to get started!**

MR CONDITIONAL/CT ELECTRODES

Rhythmink's Disposable MR Conditional/CT Electrodes are FDA Cleared for 0.2, 0.5, 1.0, 1.5 and 3T environments. Designed for more efficient work flow, these electrodes can remain on a patient safely within certain MRI environments when disconnected from their extension cables, or produce quality recordings and an artifact free image when used during CT scans.

Disposable Quick Connect System™

Designed for one box per patient use.

Product Number	Description	Lead Length	Leads Per Box						
 Slim Cups™									
MRQS2070615	2 Arrays (one 7 lead array & one 6 lead array) <table><tr><td>Array 1 (7)</td><td>Array 2 (6)</td></tr><tr><td> </td><td> </td></tr></table>	Array 1 (7)	Array 2 (6)			1.5m	13		
Array 1 (7)	Array 2 (6)								
MRQS309090515	3 Arrays (two 9 lead arrays & one 5 lead array) <table><tr><td>Array 1 (9)</td><td>Array 2 (9)</td><td>Array 3 (5)</td></tr><tr><td> </td><td> </td><td> </td></tr></table>	Array 1 (9)	Array 2 (9)	Array 3 (5)				1.5m	23
Array 1 (9)	Array 2 (9)	Array 3 (5)							
MRQS310100315	3 Arrays (two 10 lead arrays & one 3 lead array) <table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (3)</td></tr><tr><td> </td><td> </td><td> </td></tr></table>	Array 1 (10)	Array 2 (10)	Array 3 (3)				1.5m	23
Array 1 (10)	Array 2 (10)	Array 3 (3)							
MRQS32410	3 Arrays (each contains 8 leads) <table><tr><td>Array 1 (8)</td><td>Array 2 (8)</td><td>Array 3 (8)</td></tr><tr><td> </td><td> </td><td> </td></tr></table>	Array 1 (8)	Array 2 (8)	Array 3 (8)				1.0m	24
Array 1 (8)	Array 2 (8)	Array 3 (8)							
MRQS32415	3 Arrays (each contains 8 leads) <table><tr><td>Array 1 (8)</td><td>Array 2 (8)</td><td>Array 3 (8)</td></tr><tr><td> </td><td> </td><td> </td></tr></table>	Array 1 (8)	Array 2 (8)	Array 3 (8)				1.5m	24
Array 1 (8)	Array 2 (8)	Array 3 (8)							
MRQS308080815	3 Arrays (each contains 8 leads) <table><tr><td>Array 1 (8)</td><td>Array 2 (8)</td><td>Array 3 (8)</td></tr><tr><td> </td><td> </td><td> </td></tr></table>	Array 1 (8)	Array 2 (8)	Array 3 (8)				1.5m	24
Array 1 (8)	Array 2 (8)	Array 3 (8)							
MRQS32615	3 Arrays (two 8 lead arrays & one 10 lead array) <table><tr><td>Array 1 (8)</td><td>Array 2 (8)</td><td>Array 3 (10)</td></tr><tr><td> </td><td> </td><td> </td></tr></table>	Array 1 (8)	Array 2 (8)	Array 3 (10)				1.5m	26
Array 1 (8)	Array 2 (8)	Array 3 (10)							

MR CONDITIONAL/CT ELECTRODES

Disposable Quick Connect System™

Designed for one box per patient use.

Product Number	Description	Lead Length	Leads Per Box					
 Slim Cups™								
MRQS42610	4 Arrays (two 7 lead arrays & two 6 lead arrays)	1.0m	26					
	<table><tr><td>Array 1 (7)</td><td>Array 2 (7)</td><td>Array 3 (6)</td><td>Array 4 (6)</td></tr><tr><td></td><td></td><td></td><td></td></tr></table>			Array 1 (7)	Array 2 (7)	Array 3 (6)	Array 4 (6)	
Array 1 (7)	Array 2 (7)	Array 3 (6)	Array 4 (6)					
								
 Deep Cups™								
MRQD10810	1 Array of 8 leads, 10-20 labeling <table><tr><td>Array 1 (8)</td></tr><tr><td></td></tr></table>	Array 1 (8)		1.0m	8			
Array 1 (8)								
								
MRQD2060710	2 Arrays (one 6 lead array & one 7 lead array), 10-20 labeling <table><tr><td>Array 1 (6)</td><td>Array 2 (7)</td></tr><tr><td></td><td></td></tr></table>	Array 1 (6)	Array 2 (7)			1.0m	13	
Array 1 (6)	Array 2 (7)							
								
MRQD2060715	2 Arrays (one 6 lead array & one 7 lead array) <table><tr><td>Array 1 (6)</td><td>Array 2 (7)</td></tr><tr><td></td><td></td></tr></table>	Array 1 (6)	Array 2 (7)			1.5m	13	
Array 1 (6)	Array 2 (7)							
								
MRQD2070615	2 Arrays (one 7 lead array & one 6 lead array), 10-20 labeling <table><tr><td>Array 1 (7)</td><td>Array 2 (6)</td></tr><tr><td></td><td></td></tr></table>	Array 1 (7)	Array 2 (6)			1.5m	13	
Array 1 (7)	Array 2 (6)							
								
MRQD20706152	2 Arrays (one 7 lead array & one 6 lead array), 10-20 labeling <table><tr><td>Array 1 (7)</td><td>Array 2 (6)</td></tr><tr><td></td><td></td></tr></table>	Array 1 (7)	Array 2 (6)			1.5m	13	
Array 1 (7)	Array 2 (6)							
								
MRQD2070715	2 Arrays (each contain 7 leads) <table><tr><td>Array 1 (7)</td><td>Array 2 (7)</td></tr><tr><td></td><td></td></tr></table>	Array 1 (7)	Array 2 (7)			1.5m	14	
Array 1 (7)	Array 2 (7)							
								
MRQD2070810	2 Arrays (one 7 lead array & one 8 lead array), 10-20 labeling <table><tr><td>Array 1 (7)</td><td>Array 2 (8)</td></tr><tr><td></td><td></td></tr></table>	Array 1 (7)	Array 2 (8)			1.0m	15	
Array 1 (7)	Array 2 (8)							
								

MR CONDITIONAL/CT ELECTRODES

Disposable Quick Connect System™

Designed for one box per patient use.

Product Number	Description	Lead Length	Leads Per Box						
 Deep Cups™									
MRQD2090910	2 Arrays (each contains 9 leads) <table><tr><th>Array 1 (9)</th><th>Array 2 (9)</th></tr><tr><td> </td><td> </td></tr></table>	Array 1 (9)	Array 2 (9)			1.0m	18		
Array 1 (9)	Array 2 (9)								
MRQD20909102	2 Arrays (each contains 9 leads), 10-20 labeling <table><tr><th>Array 1 (9)</th><th>Array 2 (9)</th></tr><tr><td> </td><td> </td></tr></table>	Array 1 (9)	Array 2 (9)			1.0m	18		
Array 1 (9)	Array 2 (9)								
MRQD21915	2 Arrays (one 9 lead array & one 8 lead array) + 2 Single MR Deep Cups <table><tr><th>Array 1 (9)</th><th>Array 2 (8)</th></tr><tr><td> </td><td> </td></tr></table>	Array 1 (9)	Array 2 (8)			1.5m	19		
Array 1 (9)	Array 2 (8)								
MRQD2101015	2 Arrays (each contains 10 leads) <table><tr><th>Array 1 (10)</th><th>Array 2 (10)</th></tr><tr><td> </td><td> </td></tr></table>	Array 1 (10)	Array 2 (10)			1.5m	20		
Array 1 (10)	Array 2 (10)								
MRQD22015	2 Arrays (each contains 10 leads) <table><tr><th>Array 1 (10)</th><th>Array 2 (10)</th></tr><tr><td> </td><td> </td></tr></table>	Array 1 (10)	Array 2 (10)			1.5m	20		
Array 1 (10)	Array 2 (10)								
MRQD307070710	3 Arrays (each contains 7 leads), 10-20 labeling <table><tr><th>Array 1 (7)</th><th>Array 2 (7)</th><th>Array 3 (7)</th></tr><tr><td> </td><td> </td><td> </td></tr></table>	Array 1 (7)	Array 2 (7)	Array 3 (7)				1.0m	21
Array 1 (7)	Array 2 (7)	Array 3 (7)							
MRQD307080815	3 Arrays (two 8 lead arrays & one 7 lead array) <table><tr><th>Array 1 (7)</th><th>Array 2 (8)</th><th>Array 3 (8)</th></tr><tr><td> </td><td> </td><td> </td></tr></table>	Array 1 (7)	Array 2 (8)	Array 3 (8)				1.5m	23
Array 1 (7)	Array 2 (8)	Array 3 (8)							
MRQD310100315	3 Arrays (two 10 lead arrays & one 3 lead array) <table><tr><th>Array 1 (10)</th><th>Array 2 (10)</th><th>Array 3 (3)</th></tr><tr><td> </td><td> </td><td> </td></tr></table>	Array 1 (10)	Array 2 (10)	Array 3 (3)				1.5m	23
Array 1 (10)	Array 2 (10)	Array 3 (3)							

MR CONDITIONAL/CT ELECTRODES

Disposable Quick Connect System™




Designed for one box per patient use.

Product Number	Description	Lead Length	Leads Per Box								
 Deep Cups™											
MRQD32315	3 Arrays (two 9 lead arrays & one 5 lead array) <table><tr><td>Array 1 (9)</td><td>Array 2 (9)</td><td>Array 3 (5)</td></tr><tr><td> </td><td> </td><td> </td></tr></table>	Array 1 (9)	Array 2 (9)	Array 3 (5)				1.5m	23		
Array 1 (9)	Array 2 (9)	Array 3 (5)									
MRQD323154	3 Arrays (two 9 lead arrays & one 5 lead array) in green/red/white sleeves <table><tr><td>Array 1 (9)</td><td>Array 2 (9)</td><td>Array 3 (5)</td></tr><tr><td> </td><td> </td><td> </td></tr></table>	Array 1 (9)	Array 2 (9)	Array 3 (5)				1.5m	23		
Array 1 (9)	Array 2 (9)	Array 3 (5)									
MRQD40606060515	4 Arrays (three 6 lead arrays & one 5 lead array), 10-20 labeling <table><tr><td>Array 1 (6)</td><td>Array 2 (6)</td><td>Array 3 (6)</td><td>Array 4 (5)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table>	Array 1 (6)	Array 2 (6)	Array 3 (6)	Array 4 (5)					1.5m	23
Array 1 (6)	Array 2 (6)	Array 3 (6)	Array 4 (5)								
MRQD308070815	3 Arrays (two 8 lead arrays & one 7 lead array) <table><tr><td>Array 1 (8)</td><td>Array 2 (7)</td><td>Array 3 (8)</td></tr><tr><td> </td><td> </td><td> </td></tr></table>	Array 1 (8)	Array 2 (7)	Array 3 (8)				1.5m	23		
Array 1 (8)	Array 2 (7)	Array 3 (8)									
MRQD308080715	3 Arrays (two 8 lead arrays & one 7 lead array) <table><tr><td>Array 1 (8)</td><td>Array 2 (8)</td><td>Array 3 (7)</td></tr><tr><td> </td><td> </td><td> </td></tr></table>	Array 1 (8)	Array 2 (8)	Array 3 (7)				1.5m	23		
Array 1 (8)	Array 2 (8)	Array 3 (7)									
MRQD305090915	3 Arrays (two 9 lead arrays & one 5 lead array) <table><tr><td>Array 1 (5)</td><td>Array 2 (9)</td><td>Array 3 (9)</td></tr><tr><td> </td><td> </td><td> </td></tr></table>	Array 1 (5)	Array 2 (9)	Array 3 (9)				1.5m	23		
Array 1 (5)	Array 2 (9)	Array 3 (9)									
MRQD309090515	3 Arrays (two 9 lead arrays & one 5 lead array) in blue/red/black sleeves <table><tr><td>Array 1 (9)</td><td>Array 2 (9)</td><td>Array 3 (5)</td></tr><tr><td> </td><td> </td><td> </td></tr></table>	Array 1 (9)	Array 2 (9)	Array 3 (5)				1.5m	23		
Array 1 (9)	Array 2 (9)	Array 3 (5)									
MRQD32410	3 Arrays (each contains 8 leads) <table><tr><td>Array 1 (8)</td><td>Array 2 (8)</td><td>Array 3 (8)</td></tr><tr><td> </td><td> </td><td> </td></tr></table>	Array 1 (8)	Array 2 (8)	Array 3 (8)				1.0m	24		
Array 1 (8)	Array 2 (8)	Array 3 (8)									

MR CONDITIONAL/CT ELECTRODES

Disposable Quick Connect System™

Designed for one box per patient use.

Product Number	Description	Lead Length	Leads Per Box						
 Deep Cups™									
MRQD32415	3 Arrays (each contains 8 leads) <table><tr><th>Array 1 (8)</th><th>Array 2 (8)</th><th>Array 3 (8)</th></tr><tr><td></td><td></td><td></td></tr></table>	Array 1 (8)	Array 2 (8)	Array 3 (8)				1.5m	24
Array 1 (8)	Array 2 (8)	Array 3 (8)							
									
MRQD309060910	3 Arrays (two 9 lead arrays & one 6 lead array), 10-20 labeling <table><tr><th>Array 1 (9)</th><th>Array 2 (6)</th><th>Array 3 (9)</th></tr><tr><td></td><td></td><td></td></tr></table>	Array 1 (9)	Array 2 (6)	Array 3 (9)				1.0m	24
Array 1 (9)	Array 2 (6)	Array 3 (9)							
									
MRQD324102	3 Arrays (two 9 lead arrays & one 6 lead array) <table><tr><th>Array 1 (9)</th><th>Array 2 (6)</th><th>Array 3 (9)</th></tr><tr><td></td><td></td><td></td></tr></table>	Array 1 (9)	Array 2 (6)	Array 3 (9)				1.0m	24
Array 1 (9)	Array 2 (6)	Array 3 (9)							
									
MRQD310041015	3 Arrays (two 10 lead arrays & one 4 lead array) <table><tr><th>Array 1 (10)</th><th>Array 2 (4)</th><th>Array 3 (10)</th></tr><tr><td></td><td></td><td></td></tr></table>	Array 1 (10)	Array 2 (4)	Array 3 (10)				1.5m	24
Array 1 (10)	Array 2 (4)	Array 3 (10)							
									
MRQD325152	3 Arrays (two 9 lead array, one 7 lead array), 10/20 labeling <table><tr><th>Array 1 (9)</th><th>Array 2 (7)</th><th>Array 3 (9)</th></tr><tr><td></td><td></td><td></td></tr></table>	Array 1 (9)	Array 2 (7)	Array 3 (9)				1.5m	25
Array 1 (9)	Array 2 (7)	Array 3 (9)							
									
MRQD308071015	3 Arrays (one 8 lead array, one 7 lead array & one 10 lead array), 10/20 labeling <table><tr><th>Array 1 (8)</th><th>Array 2 (7)</th><th>Array 3 (10)</th></tr><tr><td></td><td></td><td></td></tr></table>	Array 1 (8)	Array 2 (7)	Array 3 (10)				1.5m	25
Array 1 (8)	Array 2 (7)	Array 3 (10)							
									
MRQD310051015	3 Arrays (two 10 lead arrays and one 5 lead array) <table><tr><th>Array 1 (10)</th><th>Array 2 (5)</th><th>Array 3 (10)</th></tr><tr><td></td><td></td><td></td></tr></table>	Array 1 (10)	Array 2 (5)	Array 3 (10)				1.5m	25
Array 1 (10)	Array 2 (5)	Array 3 (10)							
									

MR CONDITIONAL/CT ELECTRODES

Disposable Quick Connect System™

Designed for one box per patient use.

Product Number	Description	Lead Length	Leads Per Box								
 Deep Cups™											
MRQD3100510152	3 Arrays (two 10 lead arrays and one 5 lead array), 10/20 labeling <table><tr><th>Array 1 (10)</th><th>Array 2 (5)</th><th>Array 3 (10)</th></tr><tr><td> </td><td> </td><td> </td></tr></table>	Array 1 (10)	Array 2 (5)	Array 3 (10)				1.5m	25		
Array 1 (10)	Array 2 (5)	Array 3 (10)									
MRQD40707060515	4 Arrays (two 7 lead arrays, one 6 lead array & one 5 lead array), 10/20 labeling <table><tr><th>Array 1 (7)</th><th>Array 2 (7)</th><th>Array 3 (6)</th><th>Array 4 (5)</th></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table>	Array 1 (7)	Array 2 (7)	Array 3 (6)	Array 4 (5)					1.5m	25
Array 1 (7)	Array 2 (7)	Array 3 (6)	Array 4 (5)								
MRQD40606060715	4 Arrays (three arrays of 6 & one 7 lead array) <table><tr><th>Array 1 (6)</th><th>Array 2 (6)</th><th>Array 3 (6)</th><th>Array 4 (7)</th></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table>	Array 1 (6)	Array 2 (6)	Array 3 (6)	Array 4 (7)					1.5m	25
Array 1 (6)	Array 2 (6)	Array 3 (6)	Array 4 (7)								
MRQD310100515	3 Arrays (two 10 lead arrays and one 5 lead array) <table><tr><th>Array 1 (10)</th><th>Array 2 (10)</th><th>Array 3 (5)</th></tr><tr><td> </td><td> </td><td> </td></tr></table>	Array 1 (10)	Array 2 (10)	Array 3 (5)				1.0m	25		
Array 1 (10)	Array 2 (10)	Array 3 (5)									
MRQD32610	3 Arrays (two 8 lead arrays & one 10 lead array) <table><tr><th>Array 1 (8)</th><th>Array 2 (8)</th><th>Array 3 (10)</th></tr><tr><td> </td><td> </td><td> </td></tr></table>	Array 1 (8)	Array 2 (8)	Array 3 (10)				1.0m	26		
Array 1 (8)	Array 2 (8)	Array 3 (10)									
MRQD32615	3 Arrays (two 8 lead arrays & one 10 lead array) <table><tr><th>Array 1 (8)</th><th>Array 2 (8)</th><th>Array 3 (10)</th></tr><tr><td> </td><td> </td><td> </td></tr></table>	Array 1 (8)	Array 2 (8)	Array 3 (10)				1.5m	26		
Array 1 (8)	Array 2 (8)	Array 3 (10)									
MRQD308081015	3 Arrays (two 8 lead arrays & one 10 lead array), 10/20 labeling <table><tr><th>Array 1 (8)</th><th>Array 2 (8)</th><th>Array 3 (10)</th></tr><tr><td> </td><td> </td><td> </td></tr></table>	Array 1 (8)	Array 2 (8)	Array 3 (10)				1.5m	26		
Array 1 (8)	Array 2 (8)	Array 3 (10)									

MR CONDITIONAL/CT ELECTRODES

Disposable Quick Connect System™

Designed for one box per patient use.

Product Number	Description	Lead Length	Leads Per Box								
 Deep Cups™											
MRQD310100615	3 Arrays (two 10 lead arrays & one 6 lead array), 10/20 labeling <table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (6)</td></tr><tr><td> </td><td> </td><td> </td></tr></table>	Array 1 (10)	Array 2 (10)	Array 3 (6)				1.5m	26		
Array 1 (10)	Array 2 (10)	Array 3 (6)									
MRQD100610152	3 Arrays (two 10 lead arrays & one 6 lead array) <table><tr><td>Array 1 (10)</td><td>Array 2 (6)</td><td>Array 3 (10)</td></tr><tr><td> </td><td> </td><td> </td></tr></table>	Array 1 (10)	Array 2 (6)	Array 3 (10)				1.5m	26		
Array 1 (10)	Array 2 (6)	Array 3 (10)									
MRQD308100815	3 Arrays (two 8 lead arrays & one 10 lead array), 10/20 labeling <table><tr><td>Array 1 (8)</td><td>Array 2 (10)</td><td>Array 3 (8)</td></tr><tr><td> </td><td> </td><td> </td></tr></table>	Array 1 (8)	Array 2 (10)	Array 3 (8)				1.5m	26		
Array 1 (8)	Array 2 (10)	Array 3 (8)									
MRQD3081008152	3 Arrays (two 8 lead arrays & one 10 lead array), 10/20 labeling <table><tr><td>Array 1 (8)</td><td>Array 2 (10)</td><td>Array 3 (8)</td></tr><tr><td> </td><td> </td><td> </td></tr></table>	Array 1 (8)	Array 2 (10)	Array 3 (8)				1.5m	26		
Array 1 (8)	Array 2 (10)	Array 3 (8)									
MRQD309080915	3 Arrays (two 9 lead arrays & one 8 lead array), 10/20 labeling <table><tr><td>Array 1 (9)</td><td>Array 2 (10)</td><td>Array 3 (9)</td></tr><tr><td> </td><td> </td><td> </td></tr></table>	Array 1 (9)	Array 2 (10)	Array 3 (9)				1.5m	26		
Array 1 (9)	Array 2 (10)	Array 3 (9)									
MRQD310061015	3 Arrays (two 10 lead arrays & one 6 lead array), 10/20 labeling <table><tr><td>Array 1 (10)</td><td>Array 2 (6)</td><td>Array 3 (10)</td></tr><tr><td> </td><td> </td><td> </td></tr></table>	Array 1 (10)	Array 2 (6)	Array 3 (10)				1.5m	26		
Array 1 (10)	Array 2 (6)	Array 3 (10)									
MRQD42610	4 Arrays (two 7 lead arrays & two 6 lead arrays) <table><tr><td>Array 1 (7)</td><td>Array 2 (7)</td><td>Array 3 (6)</td><td>Array 4 (6)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table>	Array 1 (7)	Array 2 (7)	Array 3 (6)	Array 4 (6)					1.0m	26
Array 1 (7)	Array 2 (7)	Array 3 (6)	Array 4 (6)								

MR CONDITIONAL/CT ELECTRODES

Disposable Quick Connect System™

Designed for one box per patient use.

Product Number	Description	Lead Length	Leads Per Box			
 Deep Cups™						
MRQD42615	4 Arrays (two 7 lead arrays & two 6 lead arrays)	1.5m	26			
	Array 1 (7)			Array 2 (7)	Array 3 (6)	Array 4 (6)
MRQD40707060615	4 Arrays (two 7 lead arrays & two 6 lead arrays), 10/20 labeling	1.5m	26			
	Array 1 (7)			Array 2 (7)	Array 3 (6)	Array 4 (6)
MRQD2715	3 Arrays (two 8 lead arrays & one 9 lead array) + 2 single MR Conditional Deep Cups, 10/20 Labeling	1.5	27			
	Array 1 (8)			Array 2 (8)	Array 3 (9)	
MRQD310100715	3 Arrays (two 10 lead arrays & one 7 lead array), 10/20 labeling	1.5m	27			
	Array 1 (10)			Array 2 (10)	Array 3 (7)	
MRQD40706070715	4 Arrays (three 7 lead arrays and one 6 lead array), 10/20 labeling	1.5m	27			
	Array 1 (7)			Array 2 (6)	Array 3 (7)	Array 4 (7)
MRQD32915	3 Arrays (two 10 lead arrays & one 7 lead array) + 2 single MR Conditional Deep Cups	1.5m	27 2			
	Array 1 (10)			Array 2 (10)	Array 3 (7)	
MRQD308080915	3 Arrays (two 8 lead arrays & one 9 lead array) + 4 single MR Conditional Deep Cups	1.5m	25 4			
	Array 1 (8)			Array 2 (8)	Array 3 (9)	

MR CONDITIONAL/CT ELECTRODES

Disposable Quick Connect System™

Designed for one box per patient use.

Product Number	Description	Lead Length	Leads Per Box																																																																
 Deep Cups™																																																																			
MRQD40607060610	4 Arrays (three 6 lead arrays & one 7 lead array) + 4 single MR Conditional Deep Cup	1.0m	29																																																																
	<table><tr><td>Array 1 (6)</td><td>Array 2 (7)</td><td>Array 3 (6)</td><td>Array 4 (6)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table> <tr><td rowspan="2">MRQD329154</td><td>3 Arrays (two 10 lead arrays & one 7 lead array) + 2 MR Sticky Pads™ in Bundle</td><td rowspan="2">2.5m</td><td rowspan="2">29</td></tr> <tr><td><table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (7)</td></tr><tr><td> </td><td> </td><td> </td></tr></table><tr><td rowspan="2">MRQD310101010</td><td>3 Arrays (each contain 10 leads), 10-10 labeling</td><td rowspan="2">1.0m</td><td rowspan="2">30</td></tr><tr><td><table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (10)</td></tr><tr><td> </td><td> </td><td> </td></tr></table><tr><td rowspan="2">MRQD310101015</td><td>3 Arrays (each contain 10 leads), 10/20 labeling</td><td rowspan="2">1.5m</td><td rowspan="2">30</td></tr><tr><td><table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (10)</td></tr><tr><td> </td><td> </td><td> </td></tr></table><tr><td rowspan="2">MRQD41010050815</td><td>4 Arrays (two 10 lead arrays, one 5 lead array & one 8 lead array)</td><td rowspan="2">1.5m</td><td rowspan="2">33</td></tr><tr><td><table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (5)</td><td>Array 4 (8)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table><tr><td rowspan="2">MRQD40710071015</td><td>4 Arrays (two 10 lead arrays & two 7 lead arrays)</td><td rowspan="2">1.5m</td><td rowspan="2">34</td></tr><tr><td><table><tr><td>Array 1 (7)</td><td>Array 2 (10)</td><td>Array 3 (7)</td><td>Array 4 (10)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table></td></tr></td></tr></td></tr></td></tr></td></tr>			Array 1 (6)	Array 2 (7)	Array 3 (6)	Array 4 (6)					MRQD329154	3 Arrays (two 10 lead arrays & one 7 lead array) + 2 MR Sticky Pads™ in Bundle	2.5m	29	<table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (7)</td></tr><tr><td> </td><td> </td><td> </td></tr></table> <tr><td rowspan="2">MRQD310101010</td><td>3 Arrays (each contain 10 leads), 10-10 labeling</td><td rowspan="2">1.0m</td><td rowspan="2">30</td></tr> <tr><td><table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (10)</td></tr><tr><td> </td><td> </td><td> </td></tr></table><tr><td rowspan="2">MRQD310101015</td><td>3 Arrays (each contain 10 leads), 10/20 labeling</td><td rowspan="2">1.5m</td><td rowspan="2">30</td></tr><tr><td><table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (10)</td></tr><tr><td> </td><td> </td><td> </td></tr></table><tr><td rowspan="2">MRQD41010050815</td><td>4 Arrays (two 10 lead arrays, one 5 lead array & one 8 lead array)</td><td rowspan="2">1.5m</td><td rowspan="2">33</td></tr><tr><td><table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (5)</td><td>Array 4 (8)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table><tr><td rowspan="2">MRQD40710071015</td><td>4 Arrays (two 10 lead arrays & two 7 lead arrays)</td><td rowspan="2">1.5m</td><td rowspan="2">34</td></tr><tr><td><table><tr><td>Array 1 (7)</td><td>Array 2 (10)</td><td>Array 3 (7)</td><td>Array 4 (10)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table></td></tr></td></tr></td></tr></td></tr>	Array 1 (10)	Array 2 (10)	Array 3 (7)				MRQD310101010	3 Arrays (each contain 10 leads), 10-10 labeling	1.0m	30	<table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (10)</td></tr><tr><td> </td><td> </td><td> </td></tr></table> <tr><td rowspan="2">MRQD310101015</td><td>3 Arrays (each contain 10 leads), 10/20 labeling</td><td rowspan="2">1.5m</td><td rowspan="2">30</td></tr> <tr><td><table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (10)</td></tr><tr><td> </td><td> </td><td> </td></tr></table><tr><td rowspan="2">MRQD41010050815</td><td>4 Arrays (two 10 lead arrays, one 5 lead array & one 8 lead array)</td><td rowspan="2">1.5m</td><td rowspan="2">33</td></tr><tr><td><table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (5)</td><td>Array 4 (8)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table><tr><td rowspan="2">MRQD40710071015</td><td>4 Arrays (two 10 lead arrays & two 7 lead arrays)</td><td rowspan="2">1.5m</td><td rowspan="2">34</td></tr><tr><td><table><tr><td>Array 1 (7)</td><td>Array 2 (10)</td><td>Array 3 (7)</td><td>Array 4 (10)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table></td></tr></td></tr></td></tr>	Array 1 (10)	Array 2 (10)	Array 3 (10)				MRQD310101015	3 Arrays (each contain 10 leads), 10/20 labeling	1.5m	30	<table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (10)</td></tr><tr><td> </td><td> </td><td> </td></tr></table> <tr><td rowspan="2">MRQD41010050815</td><td>4 Arrays (two 10 lead arrays, one 5 lead array & one 8 lead array)</td><td rowspan="2">1.5m</td><td rowspan="2">33</td></tr> <tr><td><table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (5)</td><td>Array 4 (8)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table><tr><td rowspan="2">MRQD40710071015</td><td>4 Arrays (two 10 lead arrays & two 7 lead arrays)</td><td rowspan="2">1.5m</td><td rowspan="2">34</td></tr><tr><td><table><tr><td>Array 1 (7)</td><td>Array 2 (10)</td><td>Array 3 (7)</td><td>Array 4 (10)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table></td></tr></td></tr>	Array 1 (10)	Array 2 (10)	Array 3 (10)				MRQD41010050815	4 Arrays (two 10 lead arrays, one 5 lead array & one 8 lead array)	1.5m	33	<table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (5)</td><td>Array 4 (8)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table> <tr><td rowspan="2">MRQD40710071015</td><td>4 Arrays (two 10 lead arrays & two 7 lead arrays)</td><td rowspan="2">1.5m</td><td rowspan="2">34</td></tr> <tr><td><table><tr><td>Array 1 (7)</td><td>Array 2 (10)</td><td>Array 3 (7)</td><td>Array 4 (10)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table></td></tr>	Array 1 (10)	Array 2 (10)	Array 3 (5)	Array 4 (8)					MRQD40710071015	4 Arrays (two 10 lead arrays & two 7 lead arrays)	1.5m	34	<table><tr><td>Array 1 (7)</td><td>Array 2 (10)</td><td>Array 3 (7)</td><td>Array 4 (10)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table>	Array 1 (7)	Array 2 (10)	Array 3 (7)	Array 4 (10)	
Array 1 (6)	Array 2 (7)	Array 3 (6)	Array 4 (6)																																																																
MRQD329154	3 Arrays (two 10 lead arrays & one 7 lead array) + 2 MR Sticky Pads™ in Bundle	2.5m	29																																																																
	<table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (7)</td></tr><tr><td> </td><td> </td><td> </td></tr></table> <tr><td rowspan="2">MRQD310101010</td><td>3 Arrays (each contain 10 leads), 10-10 labeling</td><td rowspan="2">1.0m</td><td rowspan="2">30</td></tr> <tr><td><table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (10)</td></tr><tr><td> </td><td> </td><td> </td></tr></table><tr><td rowspan="2">MRQD310101015</td><td>3 Arrays (each contain 10 leads), 10/20 labeling</td><td rowspan="2">1.5m</td><td rowspan="2">30</td></tr><tr><td><table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (10)</td></tr><tr><td> </td><td> </td><td> </td></tr></table><tr><td rowspan="2">MRQD41010050815</td><td>4 Arrays (two 10 lead arrays, one 5 lead array & one 8 lead array)</td><td rowspan="2">1.5m</td><td rowspan="2">33</td></tr><tr><td><table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (5)</td><td>Array 4 (8)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table><tr><td rowspan="2">MRQD40710071015</td><td>4 Arrays (two 10 lead arrays & two 7 lead arrays)</td><td rowspan="2">1.5m</td><td rowspan="2">34</td></tr><tr><td><table><tr><td>Array 1 (7)</td><td>Array 2 (10)</td><td>Array 3 (7)</td><td>Array 4 (10)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table></td></tr></td></tr></td></tr></td></tr>			Array 1 (10)	Array 2 (10)	Array 3 (7)				MRQD310101010	3 Arrays (each contain 10 leads), 10-10 labeling	1.0m	30	<table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (10)</td></tr><tr><td> </td><td> </td><td> </td></tr></table> <tr><td rowspan="2">MRQD310101015</td><td>3 Arrays (each contain 10 leads), 10/20 labeling</td><td rowspan="2">1.5m</td><td rowspan="2">30</td></tr> <tr><td><table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (10)</td></tr><tr><td> </td><td> </td><td> </td></tr></table><tr><td rowspan="2">MRQD41010050815</td><td>4 Arrays (two 10 lead arrays, one 5 lead array & one 8 lead array)</td><td rowspan="2">1.5m</td><td rowspan="2">33</td></tr><tr><td><table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (5)</td><td>Array 4 (8)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table><tr><td rowspan="2">MRQD40710071015</td><td>4 Arrays (two 10 lead arrays & two 7 lead arrays)</td><td rowspan="2">1.5m</td><td rowspan="2">34</td></tr><tr><td><table><tr><td>Array 1 (7)</td><td>Array 2 (10)</td><td>Array 3 (7)</td><td>Array 4 (10)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table></td></tr></td></tr></td></tr>	Array 1 (10)	Array 2 (10)	Array 3 (10)				MRQD310101015	3 Arrays (each contain 10 leads), 10/20 labeling	1.5m	30	<table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (10)</td></tr><tr><td> </td><td> </td><td> </td></tr></table> <tr><td rowspan="2">MRQD41010050815</td><td>4 Arrays (two 10 lead arrays, one 5 lead array & one 8 lead array)</td><td rowspan="2">1.5m</td><td rowspan="2">33</td></tr> <tr><td><table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (5)</td><td>Array 4 (8)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table><tr><td rowspan="2">MRQD40710071015</td><td>4 Arrays (two 10 lead arrays & two 7 lead arrays)</td><td rowspan="2">1.5m</td><td rowspan="2">34</td></tr><tr><td><table><tr><td>Array 1 (7)</td><td>Array 2 (10)</td><td>Array 3 (7)</td><td>Array 4 (10)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table></td></tr></td></tr>	Array 1 (10)	Array 2 (10)	Array 3 (10)				MRQD41010050815	4 Arrays (two 10 lead arrays, one 5 lead array & one 8 lead array)	1.5m	33	<table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (5)</td><td>Array 4 (8)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table> <tr><td rowspan="2">MRQD40710071015</td><td>4 Arrays (two 10 lead arrays & two 7 lead arrays)</td><td rowspan="2">1.5m</td><td rowspan="2">34</td></tr> <tr><td><table><tr><td>Array 1 (7)</td><td>Array 2 (10)</td><td>Array 3 (7)</td><td>Array 4 (10)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table></td></tr>	Array 1 (10)	Array 2 (10)	Array 3 (5)	Array 4 (8)					MRQD40710071015	4 Arrays (two 10 lead arrays & two 7 lead arrays)	1.5m	34	<table><tr><td>Array 1 (7)</td><td>Array 2 (10)</td><td>Array 3 (7)</td><td>Array 4 (10)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table>	Array 1 (7)	Array 2 (10)	Array 3 (7)	Array 4 (10)														
Array 1 (10)	Array 2 (10)	Array 3 (7)																																																																	
MRQD310101010	3 Arrays (each contain 10 leads), 10-10 labeling	1.0m	30																																																																
	<table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (10)</td></tr><tr><td> </td><td> </td><td> </td></tr></table> <tr><td rowspan="2">MRQD310101015</td><td>3 Arrays (each contain 10 leads), 10/20 labeling</td><td rowspan="2">1.5m</td><td rowspan="2">30</td></tr> <tr><td><table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (10)</td></tr><tr><td> </td><td> </td><td> </td></tr></table><tr><td rowspan="2">MRQD41010050815</td><td>4 Arrays (two 10 lead arrays, one 5 lead array & one 8 lead array)</td><td rowspan="2">1.5m</td><td rowspan="2">33</td></tr><tr><td><table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (5)</td><td>Array 4 (8)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table><tr><td rowspan="2">MRQD40710071015</td><td>4 Arrays (two 10 lead arrays & two 7 lead arrays)</td><td rowspan="2">1.5m</td><td rowspan="2">34</td></tr><tr><td><table><tr><td>Array 1 (7)</td><td>Array 2 (10)</td><td>Array 3 (7)</td><td>Array 4 (10)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table></td></tr></td></tr></td></tr>			Array 1 (10)	Array 2 (10)	Array 3 (10)				MRQD310101015	3 Arrays (each contain 10 leads), 10/20 labeling	1.5m	30	<table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (10)</td></tr><tr><td> </td><td> </td><td> </td></tr></table> <tr><td rowspan="2">MRQD41010050815</td><td>4 Arrays (two 10 lead arrays, one 5 lead array & one 8 lead array)</td><td rowspan="2">1.5m</td><td rowspan="2">33</td></tr> <tr><td><table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (5)</td><td>Array 4 (8)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table><tr><td rowspan="2">MRQD40710071015</td><td>4 Arrays (two 10 lead arrays & two 7 lead arrays)</td><td rowspan="2">1.5m</td><td rowspan="2">34</td></tr><tr><td><table><tr><td>Array 1 (7)</td><td>Array 2 (10)</td><td>Array 3 (7)</td><td>Array 4 (10)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table></td></tr></td></tr>	Array 1 (10)	Array 2 (10)	Array 3 (10)				MRQD41010050815	4 Arrays (two 10 lead arrays, one 5 lead array & one 8 lead array)	1.5m	33	<table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (5)</td><td>Array 4 (8)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table> <tr><td rowspan="2">MRQD40710071015</td><td>4 Arrays (two 10 lead arrays & two 7 lead arrays)</td><td rowspan="2">1.5m</td><td rowspan="2">34</td></tr> <tr><td><table><tr><td>Array 1 (7)</td><td>Array 2 (10)</td><td>Array 3 (7)</td><td>Array 4 (10)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table></td></tr>	Array 1 (10)	Array 2 (10)	Array 3 (5)	Array 4 (8)					MRQD40710071015	4 Arrays (two 10 lead arrays & two 7 lead arrays)	1.5m	34	<table><tr><td>Array 1 (7)</td><td>Array 2 (10)</td><td>Array 3 (7)</td><td>Array 4 (10)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table>	Array 1 (7)	Array 2 (10)	Array 3 (7)	Array 4 (10)																									
Array 1 (10)	Array 2 (10)	Array 3 (10)																																																																	
MRQD310101015	3 Arrays (each contain 10 leads), 10/20 labeling	1.5m	30																																																																
	<table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (10)</td></tr><tr><td> </td><td> </td><td> </td></tr></table> <tr><td rowspan="2">MRQD41010050815</td><td>4 Arrays (two 10 lead arrays, one 5 lead array & one 8 lead array)</td><td rowspan="2">1.5m</td><td rowspan="2">33</td></tr> <tr><td><table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (5)</td><td>Array 4 (8)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table><tr><td rowspan="2">MRQD40710071015</td><td>4 Arrays (two 10 lead arrays & two 7 lead arrays)</td><td rowspan="2">1.5m</td><td rowspan="2">34</td></tr><tr><td><table><tr><td>Array 1 (7)</td><td>Array 2 (10)</td><td>Array 3 (7)</td><td>Array 4 (10)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table></td></tr></td></tr>			Array 1 (10)	Array 2 (10)	Array 3 (10)				MRQD41010050815	4 Arrays (two 10 lead arrays, one 5 lead array & one 8 lead array)	1.5m	33	<table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (5)</td><td>Array 4 (8)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table> <tr><td rowspan="2">MRQD40710071015</td><td>4 Arrays (two 10 lead arrays & two 7 lead arrays)</td><td rowspan="2">1.5m</td><td rowspan="2">34</td></tr> <tr><td><table><tr><td>Array 1 (7)</td><td>Array 2 (10)</td><td>Array 3 (7)</td><td>Array 4 (10)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table></td></tr>	Array 1 (10)	Array 2 (10)	Array 3 (5)	Array 4 (8)					MRQD40710071015	4 Arrays (two 10 lead arrays & two 7 lead arrays)	1.5m	34	<table><tr><td>Array 1 (7)</td><td>Array 2 (10)</td><td>Array 3 (7)</td><td>Array 4 (10)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table>	Array 1 (7)	Array 2 (10)	Array 3 (7)	Array 4 (10)																																				
Array 1 (10)	Array 2 (10)	Array 3 (10)																																																																	
MRQD41010050815	4 Arrays (two 10 lead arrays, one 5 lead array & one 8 lead array)	1.5m	33																																																																
	<table><tr><td>Array 1 (10)</td><td>Array 2 (10)</td><td>Array 3 (5)</td><td>Array 4 (8)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table> <tr><td rowspan="2">MRQD40710071015</td><td>4 Arrays (two 10 lead arrays & two 7 lead arrays)</td><td rowspan="2">1.5m</td><td rowspan="2">34</td></tr> <tr><td><table><tr><td>Array 1 (7)</td><td>Array 2 (10)</td><td>Array 3 (7)</td><td>Array 4 (10)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table></td></tr>			Array 1 (10)	Array 2 (10)	Array 3 (5)	Array 4 (8)					MRQD40710071015	4 Arrays (two 10 lead arrays & two 7 lead arrays)	1.5m	34	<table><tr><td>Array 1 (7)</td><td>Array 2 (10)</td><td>Array 3 (7)</td><td>Array 4 (10)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table>	Array 1 (7)	Array 2 (10)	Array 3 (7)	Array 4 (10)																																															
Array 1 (10)	Array 2 (10)	Array 3 (5)	Array 4 (8)																																																																
MRQD40710071015	4 Arrays (two 10 lead arrays & two 7 lead arrays)	1.5m	34																																																																
	<table><tr><td>Array 1 (7)</td><td>Array 2 (10)</td><td>Array 3 (7)</td><td>Array 4 (10)</td></tr><tr><td> </td><td> </td><td> </td><td> </td></tr></table>			Array 1 (7)	Array 2 (10)	Array 3 (7)	Array 4 (10)																																																												
Array 1 (7)	Array 2 (10)	Array 3 (7)	Array 4 (10)																																																																

MR CONDITIONAL/CT ELECTRODES

Disposable Quick Connect System™

Product Number	Description	Lead Length	Leads Per Box								
 Webbs™											
MRQW32410	3 Arrays (each contains 8 leads) <table><tr><th>Array 1 (8)</th><th>Array 2 (8)</th><th>Array 3 (8)</th></tr><tr><td></td><td></td><td></td></tr></table>	Array 1 (8)	Array 2 (8)	Array 3 (8)				1.0m	24		
Array 1 (8)	Array 2 (8)	Array 3 (8)									
											
MRQW32415	3 Arrays (each contains 8 leads) <table><tr><th>Array 1 (8)</th><th>Array 2 (8)</th><th>Array 3 (8)</th></tr><tr><td></td><td></td><td></td></tr></table>	Array 1 (8)	Array 2 (8)	Array 3 (8)				1.5m	24		
Array 1 (8)	Array 2 (8)	Array 3 (8)									
											
MRQW3241520	3 Arrays (two 9 lead arrays & one 5 lead array) + 1 single MR webb <table><tr><th>Array 1 (9)</th><th>Array 2 (9)</th><th>Array 3 (5)</th></tr><tr><td></td><td></td><td></td></tr></table>	Array 1 (9)	Array 2 (9)	Array 3 (5)				1.5m 2.0m	23 1		
Array 1 (9)	Array 2 (9)	Array 3 (5)									
											
MRQW308080915	3 Arrays (two 8 lead arrays & one 9 lead array), 10/20 labeling <table><tr><th>Array 1 (8)</th><th>Array 2 (8)</th><th>Array 3 (9)</th></tr><tr><td></td><td></td><td></td></tr></table>	Array 1 (8)	Array 2 (8)	Array 3 (9)				1.5m	25		
Array 1 (8)	Array 2 (8)	Array 3 (9)									
											
MRQW32610	3 Arrays (two 8 lead arrays & one 10 lead array) <table><tr><th>Array 1 (8)</th><th>Array 2 (8)</th><th>Array 3 (10)</th></tr><tr><td></td><td></td><td></td></tr></table>	Array 1 (8)	Array 2 (8)	Array 3 (10)				1.0m	26		
Array 1 (8)	Array 2 (8)	Array 3 (10)									
											
MRQW32615	3 Arrays (two 8 lead arrays & one 10 lead array) <table><tr><th>Array 1 (8)</th><th>Array 2 (8)</th><th>Array 3 (10)</th></tr><tr><td></td><td></td><td></td></tr></table>	Array 1 (8)	Array 2 (8)	Array 3 (10)				1.5m	26		
Array 1 (8)	Array 2 (8)	Array 3 (10)									
											
MRQW42615	4 Arrays (two 7 lead arrays & two 6 lead arrays) <table><tr><th>Array 1 (7)</th><th>Array 2 (7)</th><th>Array 3 (6)</th><th>Array 4 (6)</th></tr><tr><td></td><td></td><td></td><td></td></tr></table>	Array 1 (7)	Array 2 (7)	Array 3 (6)	Array 4 (6)					1.5m	26
Array 1 (7)	Array 2 (7)	Array 3 (6)	Array 4 (6)								
											

MR CONDITIONAL/CT ELECTRODES

Disposable Quick Connect System™

Designed for one box per patient use.

Product Number	Description	Lead Length	Leads Per Box									
 Webbs™												
MRQW309091015	3 Arrays (two 9 lead arrays & one 10 lead array), 10/20 labeling <table><tr><th>Array 1 (9)</th><th>Array 2 (9)</th><th>Array 3 (10)</th></tr><tr><td> </td><td> </td><td> </td></tr><tr><td> </td><td> </td><td> </td></tr></table>	Array 1 (9)	Array 2 (9)	Array 3 (10)							1.5m	29
Array 1 (9)	Array 2 (9)	Array 3 (10)										
Product Number	Description	Lead Length	Leads Per Box									
 Small Webbs™												
MRQSWB2090815	2 Arrays (one 9 lead array & one 8 lead array) <table><tr><th>Array 1 (9)</th><th>Array 2 (8)</th></tr><tr><td> </td><td> </td></tr><tr><td> </td><td> </td></tr></table>	Array 1 (9)	Array 2 (8)					1.5m	17			
Array 1 (9)	Array 2 (8)											

Small Webb™

- Our same high quality flat electrode design in a smaller size
- At just 6 mm wide, it's perfect any time you need a smaller electrode
- Perfect for delicate skin or areas on the patient where skin health or skin injury is a concern





Compare: Rhythmlink
Webb Electrode
(left) to Small Webb
Electrode (right)

MR CONDITIONAL/CT ELECTRODES

Disposable Quick Connect System™

Designed for one box per patient use.

Product Number	Description	Lead Length	Leads Per Box
Combo			
MRQDW10630	1 Arrays (6 leads) Array 1 (6) 	3.0m	6
MRQWSD308050815	3 Arrays (two 8 lead arrays & one 5 lead array) Deep Cups, Slim Cups & Webbs, 10/20 labeling Array 1 (8)  Array 2 (5)  Array 3 (8) 	1.5m	21
MRQWD307070815	3 Arrays (two 7 lead arrays & one 8 lead array) Deep Cups & Webbs, 10/20 labeling Array 1 (7)  Array 2 (7)  Array 3 (8) 	1.5m	22
MRQS232415	3 Arrays (two 9 lead arrays & one 4 lead array) + 2 Single MR Conditional Sticky Pads Array 1 (9)  Array 2 (9)  Array 3 (4) 	1.5m	22 2
MRQWSD32715	3 Arrays (two 9 lead arrays and one 7 lead array) Webbs, Slim and Deep Cups + 2 Single MR Conditional Sticky Pads Array 1 (7)  Array 2 (9)  Array 3 (9) 	1.5m	27
MRQDW32815	3 Arrays (two 10 lead arrays & one 7 lead array) Deep Cups & Webbs + 1 MR Single Deep Cup Array 1 (10)  Array 2 (10)  Array 3 (7) 	1.5m	27 1
MRQD4241525	4 Arrays (two 4 lead, one 5 lead and one 3 lead arrays) + 8 Single MR Conditional Sticky Pads Array 1 (4)  Array 2 (5)  Array 3 (4)  Array 4 (3) 	2.5m	24

MR CONDITIONAL/CT ELECTRODES

Disposable Quick Connect System™

Designed for one box per patient use.

Product Number	Description	Lead Length	Leads Per Box			
Combo						
MRQD42815	4 Arrays (two 7 lead & two 5 lead arrays) + 4 Single MR Conditional Sticky Pads	1.5m cups 2.5m pads	24 4			
	Array 1 (7)			Array 2 (5)	Array 3 (7)	Array 4 (5)
MRQD33020	3 Arrays (two 10 lead and one 8 lead array) + 2 Single MR Conditional Sticky Pads	2.0m cups 2.5m pads	28 2			
	Array 1 (10)			Array 2 (8)	Array 3 (10)	
MRQS33215	3 Arrays (two 10 lead & one 8 lead array) + 4 Single MR Conditional Sticky Pads	15m	32			
	Array 1 (10)			Array 2 (10)	Array 3 (8)	


MR CONDITIONAL/CT ELECTRODES

Disposable Individual Connectors

Product Number	Description	Lead Length	Electrodes Per Box
 Slim Cups™			
MRC210	Numbers 1 - 12 Lead Wire  Connector 	1.0m	12
MRC215		1.5m	
MRC315	Numbers 13 - 24 Lead Wire  Connector 	1.5m	
MRC415	Numbers 25 - 36 Lead Wire  Connector 	1.5m	
MRC515	Numbers 37 - 48 Lead Wire  Connector 	1.5m	
 Deep Cups™			
MRDC210	Numbers 1 - 12 Lead Wire  Connector 	1.0m	12
MRDC215		1.5m	
MRDC225		2.5m	
MRDC310	Numbers 13 - 24 Lead Wire  Connector 	1.0m	
MRDC315		1.5m	
MRDC415	Numbers 25 - 36 Lead Wire  Connector 	1.5m	20
MRDC515	Numbers 37 - 48 Lead Wire  Connector 	1.5m	
MRDC2015	Mixed Numbers Lead Wire  Connector  Number 1 3 5 7 9 10 12 14 16 18 19 20 21 22 23 24 25 27 28 29	1.5m	
MRDC2215	Numbers 1 - 22 Lead Wire  Connector 	1.5m	22
MRDC2415	3 Bundles of Numbers 1-24 Bundle 1:  Bundle 2:  Bundle 3: 		24
MRDC2615	Numbers 1 - 26 Lead Wire  Connector 		26


MR CONDITIONAL/CT ELECTRODES

Disposable Individual Connectors

Product Number	Description	Lead Length	Electrodes Per Box
 Deep Cups™			
MRDC2710	Numbers 1-27 Lead Wire  Connector 	1.0m	27
MRDC2715	One bundle of (3) bundles Deep Cup + 4 Sticky Pads 	1.5m	27
MRDC2910	Numbers 1 - 29 Lead Wire  Connector 	1.0m	29
MRDC2915		1.5m	
RLSP518	MR/CT Deep Cup with 10/20 labeling in Leadwire Sleeves Large Recording Sticky Pads™	1.0m 1.5m	25 4
MRDC3015	Numbers 1 - 30 Lead Wire  Connector 	1.5m	30
 Webbs™			
MRW615	10/20 Labeling Lead Wire  Connector 	1.5m	6
MRW215	Numbers 1 - 12 Lead Wire  Connector 	1.5m	
MRW415	Numbers 25 - 36 Lead Wire 	1.5m	
MRW2110	Numbers 1 - 21 Lead Wire  Connector 	1.0m	21
RLSP516	MR Webb Electrodes  MR Deep Cup Electrodes 	1.5m	7 15
MRW2415	Single pouch with one coil of Numbers 1-12 and one coil of Numbers 13-24 with blue and red leadwire sleeves with 10/20 labeling	1.5-2.0m	24
MRW2515	Numbers 1 - 25 Lead Wire  Connector 	1.5m	25
MRW3015	Numbers 1 - 30 Lead Wire  Connector 	1.5m	30

MR CONDITIONAL/CT ELECTRODES

Disposable Individual Connectors

Product Number	Description	Lead Length	Electrodes Per Box
 Small Webbs™			
MRSWB1415		1.5m	14
 PressOn™			
MRP210	Numbers 1 - 12 Lead Wire  Connector 	1.0m	12
MRP215	Lead Wire  Connector 	1.5m	
MRP310	Numbers 13 - 24 Lead Wire  Connector 	1.0m	
MRP315	Lead Wire  Connector 	1.5m	
MRP415	Numbers 25 - 36 Lead Wire  Connector 	1.5m	
MRPB215		1.5m	
LeadLock™			
MRLL1	1 set / pouch [set = 1 yellow & 1 red]		1
MRLL24	24 sets/ pouch [set = 1 yellow & 1 red]		24



LeadLock™
Can be used to keep leadwires organized at the disconnection point for MR Conditional/CT Individual Connectors

MR CONDITIONAL/CT ELECTRODES



RhythmLink's MR Conditional Sticky Pad™ Surface Electrodes are the same high-quality surface electrodes you've come to expect from RhythmLink, now made with improved leadwire for enhanced CT Compatibility, and FDA Clearance for 1.5 and 3T MRI. Disposable MR Conditional Sticky Pads are available in two sizes and specially formulated for Stimulating, Recording or Grounding.

MR Conditional/CT Disposable Sticky Pad Surface Electrodes

Product Number	Description	Lead Length	Electrodes Per Box
Single Stimulating Sticky Pads			
MRSPSS8015	Small [15 x 20mm] 	1.5m	80
MRSPSS8025		2.5m	
MRSPSL8015	Large [20 x 25mm] 	1.5m	
MRSPSL8025		2.5m	
Single Recording Sticky Pads			
MRSPRS1410	Small [15x20mm] 	1.0m	14
MRSPRS4025	Small [15x20mm] 	2.5m	40
MRSPRL4015	Large [20 x 25mm] 	1.5m	40
MRSPRS8015	Small [15 x 20mm] 	1.5m	80
MRSPRS8025		2.5m	
MRSPRL8015	Large [20 x 25mm] 	1.5m	
MRSPRL8025		2.5m	
Ground Sticky Pads			
MRSPG4015	45 x 35mm 	1.5m	40
MRSPG4025		2.5m	
Sample Pack			
MRSPSTMRECG1525	2 of each Sticky Pad Electrodes (Ground, Stimulating, Recording) in both Small and Large	1.5m & 2.5m	18



Ground Sticky Pad Surface Electrode



ARTIFACT FREE INVISA-ELECTRODES™ FOR CT

Have you been looking for an electrode that can be left on during a CT scan that shows no artifact? RhythmLink has the answer for you. Available in Slim Cup™, Deep Cup™, Webb™ and PressOn™ versions.

Disposable Invisa-Electrodes™

Product Number	Description	Lead Length	Electrodes Per Box
 InvisaSlim™ Cup			
CTC115		1.5m	12
CTSW2715	Slim  Webb 	1.5m	25 2
 InvisaDeep™ Cup			
CTDC115		1.5m	12
CTDC2215		1.5m	22
CTDC2615		1.5m	26
CTDC261520	Bundle 1:  Bundle 2: 	1.5m	26
CTDC26151020		1.5m	26
CTDC2715		1.5m	27
CTDC2825		2.5m	28
CTDC3020		2.0m	30
CTDC3515		1.5m	35
CTDC4415	Bundle 1:  Bundle 2: 	1.5m	44
 InvisaWebb™			
CTW115		1.5m	12
CTW2715		1.5m	27
 InvisaNeedle™ - 13mm			
CTN115		1.5m	24
CTN125		2.5m	24

RHYTHMLINK EEG TEMPLATE

Perfect for critical care environments, the RhythmLink Template is a disposable elastic template that fits over your patient's head to assist in fast and accurate placement of EEG electrodes!


The template is offered in several sizes depending on product type. Designed for universal applications, the standard template offers three varying sizes to ensure you always have the best fit for your patient. Our second option, designed for the ICU, is offered in two sizes. Both varieties are labeled with common electrode placements, 16 for ICU and 21 for standard, to ensure a fast and accurate application of electrodes.

Disposable RhythmLink Template

Product Number	Size	Qty
Standard Templates		
STDTS	48-55 cm	1
STD TM	55-60 cm	1
STD TL	60-65 cm	1
ICU Templates		
ICUTSM	45-55 cm	1
ICUTML	55-65 cm	1

Disposable RhythmLink Template Electrodes

Product Number	Includes	Qty
RLITN1415	Standard Subdermal Needles	14
RLITN1515	Standard Subdermal Needles	15
CTGN115	Guardian Needle 13mm	24
RLIT131-1.0	Deep Cup EEG Electrodes Recording Sticky Pads	21 4
RLITD2910	Deep Cup EEG Electrodes Recording Sticky Pads	25 4




RHYTHMLINK BRAINHEALTH HEADSET™

RhythmLink's BrainHealth Headset is your EEG companion for ER, ICU and critical care environments.

RhythmLink's BrainHealth Headset includes an application guide and comes pre-loaded with 21 RhythmLink PressOn™ Electrodes. The headset is easy to apply, designed to be used by medical professionals already at bedside.

Disposable RhythmLink BrainHealth™ Headset

Product Number	Electrode Diameter	Lead Length	Includes	Qty
DHP102121	7mm	1.5m	21 PressOn™ Electrodes	1



PRESSON™ ELECTRODES

Place electrodes and record data quickly with the Disposable PressOn EEG Electrode. Utilize your staff that is already present and achieve reliable, real time data to help establish the best patient care plan possible.

Disposable PressOn™ Button Electrodes

Product Number	Description	Lead Length	Quantity
MRPB215 (MR/CT)		1.5m	12



STANDARD EEG ELECTRODES



Our standard EEG electrodes provide exceptional recording quality and can be used for standard CT imaging.

Disposable Standard EEG Electrodes

Product Number	Description	Lead Length	#/Box
 Slim Cups™			
DAGD102600	Standard 10 Colors 	1.0m	10
DAGD152600		1.5m	
DAGD202600		2.0m	
DAGD252600		2.5m	
DAGD302600		3.0m	
SC1420		2.0m	14
SC2315T*		1.5m	23
SC2325		2.5m	
SC2410		1.0m	24
SC2515T*		1.5m	25
SC2520		2.0m	
SC11026		1.0m	26
SC2615		1.5m	
SC3015			30

STANDARD EEG ELECTRODES

Disposable Standard EEG Electrodes

Product Number	Description	Lead Length	#/Box
 Deep Cups™			
DC110	Standard 10 Colors 	1.0m	10
DC115		1.5m	
DC115T*		1.5m	
DC120		2.0m	
DC125		2.5m	
DC11512		1.5m	12
DC1415			14
DC1520		2.0m	15
DCPSG1520		2.0m	
DC1615		1.5m	16
DC1715		1.5m	17
SPDC1810		1.0m	18
DC21151020		1.5m	21
DC2115	 with 10/20 labeling and leadwire sleeve	1.5m	21
DC2215		1.5m	22
DC2310		1.0m	23
DC2315		1.5m	23
DC23152	 1 Bundle with 10/20 labeling		23
DC2415			24
DC24153	 with 10/20 labeling	1.5m	
DC24151020	 with 10/20 labeling	1.5m	
DC241510202	 with 10/20 labeling	1.5m	
DC2510		1.0m	25
DC2515	 with 10/20 labeling	1.5m	
DC25156	Bundle 1:  Bundle 2:  Bundle 3:  with 10/20 labeling	1.5m	
DC2520		2.0m	
DC2520T*		2.0m	
DC25101020		1.0m	
DC25151020	Bundle 1:  Bundle 2:  Bundle 3:  with 10/20 labeling	1.5m	

STANDARD EEG ELECTRODES

Disposable Standard EEG Electrodes

Product Number	Description	Lead Length	#/Box
 Deep Cups™			
DC2515201020		1.5m 2.0m	23 2
DC11526		1.5m	26
DC11526T*		1.5m	
DC2615	 with 10/20 labeling	1.5m	
DC2620		2.0m	
DC2615201020		1.5m 2.0m	24 2
DC2710		1.0m	27
DC22710		1.0m	27
DC22715		1.5m	27
DC2715		1.5m	27
DC27152	 1 bundle with 10/20 labeling in leadwire sleeve	1.5m	27
DC27102015	 with 10/20 labeling	1.5m	27
DC27151020	 with 10/20 labeling	1.5m	27
DC2720		2.0m	27
DC2815	Bundle 1:  Bundle 2:  Bundle 3:  with 10/20 labeling in leadwire sleeves	1.5m	28
DC28152	 2 Bundles of 14 electrodes with 10/20 labeling	1.5m	28
DC291520		1.5m 2.0m	27 2
DC2910		1.0m	29
DC2915	 3 Bundles in leadwire sleeves	1.5m 2.0m	29
DC3015		1.5m	30
DC31151020	Bundle 1:  Bundle 2:  Bundle 3:  with 10/20 labeling in leadwire sleeves	1.5m	31
DC33151020	 with 10/20 labeling	1.5m	33


STANDARD EEG ELECTRODES

Disposable Standard EEG Electrodes

Product Number	Description	Lead Length	#/Box
 Webbs™			
WEBB102900	Standard 10 Colors 	1.0m	10
WEBB152900		1.5m	
WEBB202900		2.0m	
WEBB252900		2.5m	
WEBB1420		2.0m	14
WEBB2315		1.5m	23
WEBB22315		1.5m	
WEBB2310C	 with Easy Pull Connectors	1.0m	
WEBB2610		1.0m	26
WEBB2615		1.5m	
WEBB2710		1.0m	27
WEBB2715		1.5m	
WEBB2715T*		1.5m	
Product Number	Description	Lead Length	#/Box
 Small Webbs™			
SWEBB1215		1.5m	12
Product Number	Description	Lead Length	#/Box
Combo			
PSGDSP1530	Deep Cups™ 	1.5m	8
	Recording Sticky Pads™  (3 pouches)	3.0m	12
DAGDWEBB1215	Slim Cups™ 	1.5m	10
	Webbs 		2
DCW1915	Deep Cups 	1.5m	10
	Webbs  with 10/20 labeling		9
RLSP517	Deep Cups 	1.5m	15
	Webbs 		7
DCDECUL22251530	Deep Cups 	2.5m	10
	Recording Sticky Pads  (3 pouches)	1.5m & 3.0m	12

STANDARD EEG ELECTRODES

Disposable Standard EEG Electrodes

Product Number	Description	Lead Length	#/Box
Combo			
SCDECUL2225	Slim Cups™  Recording Sticky Pads™  (3 pouches)	2.5m 2.5m & 3.0m	10 12
DAGDWEBB2315S	Slim Cups & Webbs™, 3 Bundles, 10/20 with Custom Colors	1.5m	23
DAGDWEBB2415	Slim Cups  Webbs 	1.5m	22 2
DAGDWEBB2615	Slim Cups  Webbs 	1.5m	24 2
DCW2615	Deep Cups™  Webbs 	1.5m	18 8
DCDECUG261525	Deep Cups  Ground Sticky Pad 	1.5m 2.5m	25 1
SCPSG1520SP2530	Slim Cups  Recording Sticky Pads (2 pouches) 	2.0m 2.5m & 3.0m	15 8
DCDECUL271015	Deep Cups  Recording Sticky Pads 	1.0m 1.5m	23 4
RLSP511	Deep Cups  Recording Sticky Pads 	1.0m 1.5m	23 4
DCW2715	Deep Cups  Webbs  with 10/20 labeling	1.5m	16 11
DCDECUL2815	Deep Cups  Recording Sticky Pads 	1.5m	24 4
DAGDECUL2915	Slim Cups  Recording Sticky Pads 	1.5m 2.5m	25 4
DCDECUL3015	Deep Cups  Recording Sticky Pads 	1.5m	26 4



EMG/NCV Products

EMG NEEDLE ELECTRODES

Disposable Detachable Monopolar Needle Electrodes for EMG

Disposable sterile detachable monopolar needle electrodes for EMG feature a PTFE coated stainless steel needle. Each package includes two 251-24TP reusable 0.61 meter (24") flexible leadwires ending with a touch proof connector.

Product Number	Description	Needle Length	Electrodes Per Box
110-725	0.35mm Diameter (27g) - Green	25mm	35
110-637	0.40mm Diameter (26g) - Pink	37mm	
110-737	0.35mm Diameter (27g) - Blue	37mm	
110-744	0.35mm Diameter (27g) - Lavender	44mm	
110-650	0.40mm Diameter (26g) - White	50mm	
251-24TP	Reusable 24" Touch Proof Leadwire		1



Disposable Monopolar Needle Electrodes for EMG

Disposable sterile PTFE coated stainless steel monopolar needle electrodes for EMG have a 0.61 meter (24") flexible PVC insulated leadwire ending with a standard touch proof connector.

Product Number	Description	Needle Length	Electrodes Per Box
111-725-24TP	0.35mm Diameter (27g) - Green	25mm	50
111-637-24TP	0.40mm Diameter (26g) - Pink	37mm	
111-737-24TP	0.35mm Diameter (27g) - Blue	37mm	
111-837-24TP	0.32mm Diameter (28g) - Blue	37mm	
111-744-24TP	0.35mm Diameter (27g) - Lavender	44mm	
111-650-24TP	0.40mm Diameter (26g) - White	50mm	
111-575-24TP	0.45mm Diameter (25g) - Gray	75mm	25



Disposable Injectable Monopolar Needle Electrodes for EMG

Disposable sterile PTFE coated stainless steel injectable monopolar needle electrodes combine EMG with related drug therapy injections such as Botulinum toxin. They have an attached 0.61 meter (24") leadwire made of PVC insulated wire ending with a standard touch proof connector. Choose box of 10 or box of 25.

Product Number	Description	Needle Length	Electrodes Per Box
242-025-24TP	0.30mm Diameter (30g) - Blue	25mm	10
242-025-24TP-25	0.30mm Diameter (30g) - Blue		25
242-725-24TP	0.41mm Diameter (27g) - Green		10
242-725-24TP-25	0.41mm Diameter (27g) - Green		25
242-637-24TP	0.46mm Diameter (26g) - Pink	37mm	10
242-637-24TP-25	0.46mm Diameter (26g) - Pink		25
242-550-24TP	0.51mm Diameter (25g) - White	50mm	10
242-550-24TP-25	0.51mm Diameter (25g) - White		25
242-375-24TP	0.64mm Diameter (23g) - Gray	75mm	10
242-3100-24TP	0.64mm Diameter (23g) - Natural	100mm	10
242-3125-24TP	0.64mm Diameter (23g) - Natural	125mm	10



EMG/NCV SURFACE ELECTRODES

Disposable Detachable Concentric Needle Electrodes for EMG

Rhythmink's quality and consistent availability are now available for EMG/NCV users looking for a concentric needle!

Product Number	Description	Needle Length	Electrodes Per Box
EMGCNR25	0.30mm Diameter (30g) - Red	25mm	20
EMGCNY25	0.40mm Diameter (27g) - Yellow	25mm	
EMGCNG37	0.45mm Diameter (26g) - Dark Green	37mm	
EMGCNB50	0.45mm Diameter (26g) - Dark Blue	50mm	
EMGCNV75	0.60mm Diameter (23g) - Purple	75mm	
EMGCNCOMBO	Combo Pack that Contains 4 of Each Concentric Needle Sizes	Various	



Reusable Cables for Concentric Needle Electrodes

Rhythmink's reusable cable for our concentric needle electrodes is a perfect fit.

Product Number	Description	Lead Length	Electrodes Per Box
EMGCN-CBL15	5-pole DIN Connector	1.5m	1



Disposable NCV Sticky Pad Surface Electrodes without Leadwire

Our disposable NCV Sticky Pad Surface Electrodes are cleared for recording and available in a variety of shapes and sizes, all with adhesive hydrogel with keeps the electrodes in place and enhances recordings.

Product Number	Description	Electrodes Per Box
NCVB100	12x18mm with 25mm between Disc Areas	100
NCVR100	8x95mm Ring Electrodes	100
NCVT100	22x32mm Tab Electrodes	100
NCVG50	40x50mm Ground Electrodes	50
NCVKIT15	In each pouch: 2 Paired Rings, 2 Bars, 2 Tabs In each box: 15 pouches	90
256-2A1J-48STP	Alligator Clips (3) to Touch Proof (Female) Gray	1 Clip
256-A-24BTP	Alligator Clip - 24" Black Leadwire	
256-A-24GTP	Alligator Clip - 24" Green Leadwire	
256-A-24RTP	Alligator Clip - 24" Red Leadwire	
256-A-36BTP	Alligator Clip - 36" Black Leadwire	
256-A-36GTP	Alligator Clip - 36" Green Leadwire	
256-A-36RTP	Alligator Clip - 36" Red Leadwire	
256-A-48BTP	Alligator Clip - 48" Black Leadwire	
256-A-48GTP	Alligator Clip - 48" Green Leadwire	
256-A-48RTP	Alligator Clip - 48" Red Leadwire	

EMG/NCV SURFACE ELECTRODES

Disposable NCV Sticky Pad Surface Electrodes with Leadwire


Our disposable NCV Sticky Pad Surface Electrodes are cleared for recording and available in a variety of shapes and sizes, all with adhesive hydrogel with keeps the electrodes in place and enhances recordings.

Product Number	Description	Lead Length	Electrodes Per Box
NCVD80L	20mm Disc Electrodes	1.5m	80
NCVT80L	20x25mm Tab Electrodes	1.0m	80
NCVR40L	8x106mm Ring Electrodes	0.5m	40
NCVG40L - 1.0	40x50mm Ground Electrodes	1.0m	40
NCVG40L - 2.0	40x50mm Ground Electrodes	2.0m	40
NCVKIT60L	2 Discs and 1 Ground Electrode Kit	1.0m	60

NCV Sticky Pads without leadwire



NCVKIT60L - 2 Discs and 1 Ground Electrode Kit



Extension Cables

Single extension cables allow you to extend the leadwire length of a standard EMG electrode. It is a 0.61 meter (24") PVC insulated leadwire ending with a standard touch proof connector. Other lengths available upon request.

Paired cables allow you to extend the leadwire length of a standard EMG electrode. It is available in 0.61 meter (24") or 1.2 meter (48") PVC insulated red and black twisted leadwire.

Product Number	Description	Leadwire	Quantity
Single			
253-24BTP	Single - Black	0.61m	1
253-24RTP	Single - Red		
253-24GTP	Single - Green		
Paired			
254-24TP	Twisted - Red & Black	0.61m	1
254-48TP	Twisted - Red & Black	1.2m	





RhythmLink - IONM Products


STIMULATING PROBES

Rhythmink's Stimulating Probes are made from stainless steel with handles specifically designed for a surgeons needs. Each probe has been designed for specific procedures in mind.

Disposable Stimulating Probes

Product Number	Probe Length	Probe Diameter	Tip Diameter	Lead Length	Quantity
Standard Ball Tip Stimulating Probe					
RLBTP101-2.5	100mm	1.5mm	3mm	2.5m	1
RLBTP100-2.5					10
Tapered Ball Tip Stimulating Probe					
RLBTP151-2.5	150mm	1.85-2.5mm	2.3mm	2.5m	1
RLBTP150-2.5					10
RLBTP201-2.5	200mm				1
RLBTP200-2.5					10
Standard Monopolar Stimulating Probe					
RLMSP051-2.5	100mm	1.5mm	0.75mm	2.5m	1
RLMSP050-2.5					10
Extended Monopolar Stimulating Probe					
RLMSP101-2.5	160mm	3.5mm	1.6mm	2.5m	1
RLMSP100-2.5					5
Monopolar Stimulating Probe with Removable Handle					
RLSP400	200mm	1.5mm	1.5mm	2.5m	1

These products are not a part of Rhythmink's sample program.



RHYTHMLINK SUBDERMAL NEEDLES



Single Subdermal Needles - Quality recording and stimulating can only occur when you're using a quality product. Our single straight disposable needles are available in the following lengths:


7mm

13mm

19mm

22mm

Single bent subdermal needles are also available.



Parallel Paired Subdermal Needles - flexibility to separate the leadwire as much or as little as needed, our paired needles are available in:

7mm

13mm

19mm

22mm




Corkscrew Subdermal Needles - designed to stay securely in place during application. Our disposable corkscrew needles offer an easy to grip hub and are packaged in individual sterile pouches.

Guardian Needle™ - designed with enhanced features to protect the patient and medical personnel from accidental needle injury as well as delivering equivalent recording and stimulating quality.




Guardian+ Plus Needle™ The Guardian+ Plus Needle encloses the hub in a sterile adhesive bandage designed to replace the use of tape to secure the needle in place.


Twisted Subdermal Needles - multiple needle electrodes offered in one, convenient easy open pouch with twisted wires for better organization. Available as straight needles in twisted pairs, twisted bundles or corkscrews in pairs or 4-Twist.



Twisted Pairs
Subdermal Needle
Electrodes



5-Twist Subdermal
Needle Electrodes



6-Twist Subdermal
Needle Electrodes

SUBDERMAL NEEDLES

RhythmLink's Disposable Subdermal Needle Electrodes offer the quality you expect and deliver the patient comfort that you need. Made of high-tensile strength stainless steel and cleared for both stimulating and recording. Select from a variety of needle lengths and styles, leadwire colors and lengths, singles or pairs.

Disposable Single Needle Electrodes

Product Number	Description	Lead Length	Electrodes Per Box
Single 7mm (Short) Needles			
RLSND107-1.5	6-Color Group 1 	1.5m	24
RLSND107-2.5		2.5m	
RLSND109-1.5	6-Color Group 2 	1.5m	24
RLSND109-2.5		2.5m	
RLSND108-2.5	12-Color 	2.5m	24
Single 13mm Needles*			
RLSND314-2.5	Multipack with Custom Colors	2.5m	14
RLSP351	5-Color Group 1 	1.5m	25
RLSP450		2.5m	20
RLSP350	5-Color Group 2 	1.5m	25
RLSND110-1.0	6-Color Group 1 	1.0m	24
RLSND110-1.5		1.5m	
RLSND110-2.5		2.5m	
RLSND112-2.5	6-Color Group 2 	2.5m	24
RLSND111-1.0	12-Color 	1.0m	24
RLSND111-1.5		1.5m	
RLSND111-2.5		2.5m	
* 13mm Needles are also available in individual color boxes, in 1.5 & 2.5m length in quantities of 24 electrodes.			
Single 19mm Needles			
RLSND113-1.5	6-Color Group 1 	1.5m	24
RLSND113-2.5		2.5m	
RLSND114-1.5	12-Color 	1.5m	24
RLSND114-2.5		2.5m	
Single 22mm Needles			
RLSND170-1.5	6-Color Group 1 	1.5m	24
RLSND170-2.5		2.5m	

SUBDERMAL NEEDLES

Disposable Paired Needle Electrodes

Product Number	Description	Lead Length	Electrodes Per Box
Paired 7mm (Short) Needles			
RLSND125-1.5	Color Group 1 	1.5m	10
RLSND125-2.5		2.5m	
RLSND126-1.5	Color Group 2 	1.5m	10
RLSND126-2.5		2.5m	
Paired 13mm Needles*			
RLSND121-1.0	Color Group 1 	1.0m	10
RLSND121-1.5		1.5m	
RLSND121-2.0		2.0m	
RLSND121-2.5		2.5m	
RLSP310		1.5m	20
RLSND122-1.0	Color Group 2 	1.0m	10
RLSND122-1.5		1.5m	
RLSND122-2.5		2.5m	
RLSP503	Color Group 4 	2.5m	12
RLSP531C	 Easy pull touchproof connectors	2.5m	12
RLSND116P-1.5	Color Group 1  13mm Paired Needle Electrodes Bent Needle Electrodes	1.5m	10
* 13mm Paired Needles are also available in individual color boxes, in 1.5 & 2.5m length in quantities of 10 electrodes.			
Paired 19mm Needles			
RLSND123-1.5	Color Group 1 	1.5m	10
RLSND123-2.5		2.5m	
RLSND124-1.5	Color Group 2 	1.5m	10
RLSND124-2.5		2.5m	
Paired 22mm Needles			
RLSND173-1.5	Color Group 1 	1.5m	10
RLSND173-2.5		2.5m	
RLSND174-1.5	Color Group 2 	1.5m	10
RLSND174-2.5		2.5m	
RLSP504	Color Group 3 	2.5m	8

Disposable Twisted Paired Needle Electrodes

Product Number	Description	Lead Length	Electrodes Per Box
Twisted Paired 13mm Needles			
RLSND131-2.5	Color Group 1 ●○ ●●●●●●	2.5m	10 pairs
RLSND132-2.5	Color Group 2 ●○ ●●●●●●	2.5m	10 pairs
RLSP430	Color Group 3 ●○ ●●●●●●	1.5m	12 pairs
TSN2325		2.5m	
TSN2320	●●●●●●●●●●●●●●	2.0m	12 pairs

Disposable 5 & 6 Twisted Needle Electrodes

Product Number	Description	Lead Length	Set of Twists Per Box
RLSP484	5-Twist ●●●●●	1.5m	4
RLSP513		2.5m	
RLSP483	6-Twist ●●●●●●	1.5m	3
RLSP512		2.5m	
RLSND118-1.5	6-Twist with Bent Needles ●●●●●●	1.5m	3
RLSP527	7-Twist	2.5m	3

Disposable Bent Needle Electrodes

Product Number	Description	Lead Length	Electrodes Per Box
SNB1710	Bundle 1 ●○●●●●●● Bundle 2 ●○●●●●●●	1.0m	17
RLSND116-1.5	6-Color Group 1 ●○●●●●	1.5m	24
RLSND116-2.5		2.5m	
RLSND117-1.5	12-Color ●○●●●●●●●●●●	1.5m	24

SUBDERMAL NEEDLES

Disposable Corkscrew Needle Electrodes

Product Number	Description	Lead Length	Electrodes Per Box
RLSP486	4 Color Single ○●●●	1.5m	12
RLSND140-1.0	6 Color Single ●○●●●●	1.0m	12
RLSND140-1.5		1.5m	
RLSND140-2.5		2.5m	
TCN2115	2-Twist Corkscrew ●●	1.5m	12
TCN4115	4-Twist Corkscrew ●○●●		12
TCN2125	2-Twist Corkscrew ●●	2.5m	12
TCN4125	4-Twist Corkscrew ●○●●	2.5m	12
TCN5115	5-Twist Corkscrew ●●●●●	1.5m	15
TCN5125	5-Twist Corkscrew ●●●●●	2.5m	15

FIXED DISTANCE GUARDIAN NEEDLE

RhythmLink's Fixed Distance Guardian Needles allow you to place two needles at once, perfect for when you need reliable and precise distance between two electrodes. The Fixed Distance Guardian Needle includes two 13mm Subdermal Needles at a fixed distance of 2.5mm, packaged with additional individual needles for convenience.

Fixed Distance Guardian Needle and Guardian+ Plus Needle

Product Number	Description	Lead Length	Electrodes Per Box
GNFD5152	3 Fixed Distance Guardian Needles and 2 Single Guardian Needles	1.5m	5
GNFD515	3 Fixed Distance Guardian+ Plus Needles and 2 Single Guardian+ Plus Needles [with Bandage]	1.5m	5



RHYTHMLINK GUARDIAN NEEDLE™

The Guardian Needle Electrode is a subdermal needle electrode which can be used in any situation where the current subdermal needles are used. Enhanced features are designed to protect the patient and medical personnel from accidental needle injury as well as delivering equivalent recording and stimulating quality.

Disposable Guardian Needle - Safety Subdermal Electrode

Product Number	Description	Lead Length	Electrodes Per Box
CTGN215	13mm Guardian Needle	1.5m	24
CTGN2121-1.5	13mm Paired Guardian Needle	1.5m	10 Pairs
CTGN2121-2.5	13mm Paired Guardian Needle	2.5m	10 Pairs



Disposable Guardian+ Plus Needle - Safety Subdermal Electrode

The Guardian+ Plus Needle offers additional features by enclosing the hub in a sterile adhesive bandage designed to replace the use of tape to secure the needle in place.

Product Number	Description	Lead Length	Electrodes Per Box
CTGN615	13mm Guardian Plus Needle [with Bandage]	1.5m	6
CTGN115		1.5m	24
CTGN121-1.5	13mm Paired Guardian Plus Needle [with Bandage]	1.5m	10 Pairs
CTGN121-2.5	13mm Paired Guardian Plus Needle [with Bandage]	2.5m	10 Pairs
RLITCTGNM1.0	13mm 10/20 Guardian+ Plus Needles [with Bandage] Standard Medium RLI Template	1.0m	21 1
RLITCTGNL1.0	13mm 10/20 Guardian+ Plus Needles [with Bandage] Standard Large RLI Template	1.0m	21 1

STICKY PAD™ SURFACE ELECTRODES

RhythmLink's Sticky Pad Surface Electrodes are disposable, pre-gelled and formulated for recording, stimulation or grounding.

Disposable Standard Sticky Pad Surface Electrodes

Product Number	Description	Lead Length	Electrodes Per Box
Single Stimulating Sticky Pads			
STCUS15026	Small [15 x 20mm] ●●●●	1.5m	80
STCUS25026		2.5m	
STCUS30026		3.0m	
STCUL15026	Large [20x25mm] ●●●●	1.5m	
STCUL25026		2.5m	
Parallel Pair Stimulating Sticky Pads			
PSTCUS30026	Small [15 x 20mm] ●●●●	3.0m	40 pairs
PSTCUS15026	Small [15 x 20mm] ●●●●	1.5m	
PSTCUS25026		2.5m	
PSTCUL15026	Large [20x25mm] ●●●●	1.5m	
PSTCUL25026		2.5m	
Single Recording Sticky Pads			
DECUS10026	Small [15 x 20mm] ●●●●	1.0m	80
DECUS15026		1.5m	
DECUS25026		2.5m	
DECUL15026	Large [20x25mm] ●●●●	1.5m	
DECUL25026		2.5m	
RLSP537	Large [20x25mm] ●●●●●	3.0m	80
RLSP492	Small [15 x 20mm] ●●●●	1.5m	60
RLSP536	Large [20x25mm] ●●●●●●●	3.0m	60
Parallel Pair Recording Sticky Pads			
PDECUS15026	Small [15 x 20mm] ●●●●	1.5m	40 pairs
PDECUS25026		2.5m	
PDECUL15026	Large [20x25mm] ●●●●	1.5m	
PDECUL25026		2.5m	
RLSP482	Small [15 x 20mm] ●●	3.0m	20 pairs
Ground Sticky Pads			
DECUG10026	45 x 35mm ●	1.0m	40
DECUG15026		1.5m	
DECUG25026		2.5m	
RLSPGSP120	40x50mm ●	2.0m	
NCV115	40 Small Recording Sticky Pads ●●●● 20 Ground Sticky Pads ●	1.5m	60

ALLIGATOR CLIPS

Rhythmink's Disposable Alligator Clips are 100% disposable, sterile and intended for use in the recording of intraoperative neuromonitoring.

Disposable Alligator Clip

Product Number	Description	Lead Length	Quantity
RLSP487	Small Clip	2.5m	12




ET TUBE ELECTRODES

Designed to easily fit any ET tube, the Dragonfly® Laryngeal Electrode is designed to accurately monitor the Xth Nerve/ RLN monitoring. With unique stripes for easy placement, the Dragonfly is compatible with most recording devices. The Disposable Cobra® Monitoring ET Tube is constructed of PVC for quality performance and provides monitoring of the RLN using a standard ET Tube. The Cobra Electrode's recording surface is flush with the ET tube for smooth airway contact.

Disposable ET Tube Electrodes

Product Number	Description	ET Tube Size	Lead Length	Quantity
Disposable Dragonfly Electrodes				
LSE500MS-5	Single Channel	6.0 - 7.0	0.8m	5
LSE500M-5		7.5 - 10.0		
LSE500DCS-5	2-Channel	6.0 - 7.0	2.0m	5
LSE500DCL-5		8.0 - 9.5		
Disposable Cobra Monitoring ET Tube				
LTE7003PS-5	3-Plate EMG Endotracheal Tube	6.0	2.0m	5
LTE7003PM-5		7.0		
LTE7003PL-5		8.0		

This product is not part of Rhythmink's sample program.
Dragonfly and Cobra are registered trademarks of Neurovision Medical Products.



SIGHTSAVER™ FLASH VISUAL EVOKED POTENTIAL STIMULATOR

The SightSaver™ is a FDA approved device to facilitate the monitoring of visual pathways, and prevent vision loss during surgery. It consists of a disposable adhesive foam padding shaped to the contours of the periocular region attached to a printed circuit board with specially selected high intensity light emitting diodes (LED). The FDA approved SightSaverMini™ is a more streamlined version of the standard SightSaver™ visor. A reusable cable interface is compatible with standard clinical neurophysiology systems used for triggering an LED goggle device.

The SightSaver

Product Number	Compatibility	Use
SightSaver Visor		
ATv001	See Cables Below	Disposable
SightSaver Visor Cables		
ATcc002	Cadwell w/ID Chip	Reusable
ATcne001	Natus Eneavor	
ATcm001	Medtronic NIM Eclipse	
ATcnk002	Nihon Kohden	
ATcnp001	Natus/Xitek Protector	
ATcutp001	Universal/Touchproof	
SightSaver Mini		
ATmv001 (one pair)	See Cables Below	Disposable
SightSaver Mini Cables		
ATmvc001		Reusable

The SightSaverVisor and SightSaverMini are not part of Rhythmlink's Sample Program.
SightSaver is a trademark of Anshel Technology Inc.



GELS & PASTES

Ten20® Conductive Paste is a washable, non-drying adhesive perfect for applying Rhythmlink's EEG Electrodes. Ten20 can also be used with collodion giving you the advantage of not having to re-gel the electrode once it's in place.


Gels & Pastes

Product Number	Size	Quantity
10-20-4	4 oz	3 Jars
10-20-4T	4 oz	3 Tubes
10-20-8	8 oz	3 Jars
10-30	4 oz	3 Tubes

Ten20 is a registered trademark of Weaver and Company.

RHYTHMLINK CUSTOM ELECTRODE OPTIONS

CUSTOM PRODUCTS - Have an idea for a new product that you'd like to see? RhythmLink has an accomplished product design, development and manufacturing team that offer years of experience for new products. No matter where you are in the design process, RhythmLink can help you get there, from concept to completed, packaged product.



PRIVATE LABELING - Select from a variety of options to have our product privately labeled for your facility or brand. RhythmLink can privately label an existing product or create a new packaged product to your specifications.



BUNDLED PACKAGING - RhythmLink can fulfill your need for personalized packaging! Whether it's more or less electrodes per pouch, multiple products in a box or single color leadwire, we can develop a solution to answer your request. Take advantage of existing product options or send in your own ideas.

Ready to make your own truly custom product?
Contact sales@rhythmink.com to get started!

PRODUCT LISTING BY PART NUMBER

PART NUMBER	PAGE				
10-30	45	DAGDWEBB2615	30	MRDC210	20
10-20-4	45	DC110	27	MRDC215	20
10-20-4T	45	DC115	27	MRDC225	20
10-20-8	45	DC115T	27	MRDC310	20
110-637	32	DC11512	27	MRDC315	20
110-650	32	DC11526	28	MRDC415	20
110-725	32	DC11526T	28	MRDC515	20
110-737	32	DC1415	27	MRDC2015	20
110-744	32	DC120	27	MRDC2215	20
111-575-24TP	32	DC125	27	MRDC2415	20
111-637-24TP	32	DC1615	27	MRDC2615	20
111-650-24TP	32	DC1520	27	MRDC2710	21
111-725-24TP	32	DC1715	27	MRDC2715	21
111-737-24TP	32	DC2115	27	MRDC2910	21
111-744-24TP	32	DC2215	27	MRDC2915	21
111-837-24TP	32	DC2310	27	MRDC3015	21
242-025-24TP	32	DC2315	27	MRLL1	22
242-025-24TP-25	32	DC23152	27	MRLL24	22
242-375-24TP	32	DC2415	27	MRP210	22
242-550-24TP	32	DC24153	27	MRP215	22
242-550-24TP-25	32	DC24151020	27	MRP310	22
242-637-24TP	32	DC241510202	27	MRP315	22
242-637-24TP-25	32	DC2510	27	MRP415	22
242-725-24TP	32	DC2515	27	MRPB215	22
242-725-24TP-25	32	DC25156	27	MRSPG4015	23
242-3100-24TP	32	DC2520	27	MRSPG4025	23
242-3125-24TP	32	DC2520T	27	MRSPRS1410	23
251-24TP	32	DC25101020	27	MRSPRS4025	23
253-24BTP	34	DC25151020	27	MRSPRS8015	23
253-24GTP	34	DC2515201020	28	MRSPRS8025	23
253-24RTP	34	DC2615	28	MRSPRL4015	23
254-24TP	34	DC2620	28	MRSPRL8015	23
254-48TP	34	DC2615201020	28	MRSPRL8025	23
256-2AIJ-48STP	33	DC2710	28	MRSPSS8015	23
256-A-24BTP	33	DC2715	28	MRSPSS8025	23
256-A-24GTP	33	DC22715	28	MRSPSL8015	23
256-A-24RTP	33	DC27151020	28	MRSPSL8025	23
256-A-36BTP	33	DC2720	28	MRSPSTMREC61525	23
256-A-36GTP	33	DC2815	28	MRSWB1415	22
256-A-36RTP	33	DC28152	28	MRQD10810	8
256-A-48BTP	33	DC2910	28	MRQD2060710	8
256-A-48GTP	33	DC2915	28	MRQD2070615	8
256-A-48RTP	33	DC3015	28	MRQD2070810	8
Atv001	45	DC31151020	28	MRQD2070715	8
ATcc002	45	DC33151020	28	MRQD2090910	9
ATcne001	45	DCDECUG261525	30	MRQD2101015	9
ATcm001	45	DCDECUL271015	30	MRQD21915	9
ATcnk002	45	DCDECUL2815	30	MRQD22015	9
ATcnp001	45	DCDECUL3015	30	MRQD307070710	9
ATcutp001	45	DCW2615	30	MRQD307080815	9
ATmv001	45	DCW2715	30	MRQD32315	10
ATmvc001	45	DECUG10026	43	MRQD323154	10
CTC115	24	DECUG15026	43	MRQD40606060515	10
CTDC115	24	DECUG25026	43	MRQD32410	10
CTDC2215	24	DECUS10026	43	MRQD32415	11
CTDC2615	24	DECUS15026	43	MRQD32715	14
CTDC261520	24	DECUS25026	43	MRQD305090915	10
CTDC26151020	24	DECUL15026	43	MRQD308070815	10
CTDC2715	24	DECUL25026	43	MRQD308080715	10
CTDC2825	24	DCPSG1520	27	MRQD308071015	11
CTDC3020	24	DHP102121	25	MRQD309090515	10
CTDC3515	24	EMGCNR25	33	MRQD309060910	11
CTDC4415	24	EMGCNY25	33	MRQD324102	11
CTGN615	42	EMGCNG37	33	MRQD310041015	11
CTGN115	42	EMGCNB50	33	MRQD40707060515	12
CTGN215	42	EMGCNV75	33	MRQD310051015	11
CTGN2121-1.5	42	EMGCN-CBL15	33	MRQD100510152	12
CTGN2121-2.5	42	EMGCNCOMBO	33	MRQD310100515	12
CTGN121-1.5	42	GNFD5152	41	MRQD40706070715	14
CTGN121-2.5	42	GNFD515	41	MRQD325152	11
CTN115	24	ICUTML	25	MRQD32610	12
CTN125	24	ICUTSM	25	MRQD32615	12
CTSW2715	24	LSE500DCS-5	44	MRQD308081015	12
CTW115	24	LSE500DCL-5	44	MRQD310100615	13
CTW2715	24	LSE500M-5	44	MRQD40606060715	12
DAGD102600	26	LSE500MS-5	44	MRQD308100815	13
DAGD152600	26	LTE7003PS-5	44	MRQD3081008152	13
DAGD202600	26	LTE7003PM-5	44	MRQD309080915	13
DAGD252600	26	LTE7003PL-5	44	MRQD310061015	13
DAGD302600	26	MRC210	20	MRQD3100610152	13
DAGDDEECUL2915	30	MRC215	20	MRQD42610	13
DAGDWEBB2315S	30	MRC315	20	MRQD42615	14
DAGDWEBB2415	30	MRC415	20	MRQD40707060615	14
		MRC515	20	MRQD310100715	14

PRODUCT LISTING BY PART NUMBER

MRQD32915	14	RLSND110-2.5	38	SWEBB1215	29
MRQD329154	15	RLSND111-1.0	38	TCN2115	41
MRQD308080915	14	RLSND111-1.5	38	TCN2125	41
MRQD310101010	15	RLSND111-2.5	38	TCN4115	41
MRQD310101015	15	RLSND112-2.5	38	TCN4125	41
MRQD41010050815	15	RLSND113-1.5	38	TCN5115	41
MRQD40710071015	15	RLSND113-2.5	38	TCN5125	41
MRQD40607060610	15	RLSND114-1.5	38	TSN2320	40
MRQD33020	19	RLSND114-2.5	38	TSN2325	40
MRQD4241525	18	RLSND116-1.5	40	WEBB102900	29
MRQD42815	19	RLSND116-2.5	40	WEBB152900	29
MRQDW10630	18	RLSND116P-1.5	39	WEBB202900	29
MRQDW32815	18	RLSND117-1.5	40	WEBB252900	29
MRQS2070615	7	RLSND118-1.5	40	WEBB1420	29
MRQS323415	18	RLSND121-1.0	39	WEBB22315	29
MRQS32410	7	RLSND121-1.5	39	WEBB2315	29
MRQS32415	7	RLSND121-2.0	39	WEBB2310C	29
MRQS32615	7	RLSND121-2.5	39	WEBB2610	29
MRQS309090515	7	RLSND122-1.0	39	WEBB2615	29
MRQS310100315	7	RLSND122-1.5	39	WEBB2710	29
MRQS308080815	7	RLSND122-2.5	39	WEBB2715	29
MRQS33215	19	RLSND123-1.5	39	WEBB2715T	29
MRQS42610	8	RLSND123-2.5	39		
MRQSWB2090815	17	RLSND124-1.5	39		
MRQW32410	16	RLSND124-2.5	39		
MRQW32415	16	RLSND125-1.5	39		
MRQW3241520	16	RLSND125-2.5	39		
MRQW308080915	16	RLSND126-1.5	39		
MRQW32610	16	RLSND126-2.5	39		
MRQW32615	16	RLSND131-2.5	40		
MRQW309091015	16	RLSND132-2.5	40		
MRQW42615	16	RLSND140-1.0	41		
MRQW309091015	17	RLSND140-1.5	41		
MRQWSD308050815	18	RLSND140-2.5	41		
MRQWSD32715	18	RLSND170-1.5	38		
MRQWD307070815	18	RLSND170-2.5	38		
MRW615	21	RLSND173-1.5	39		
MRW2110	21	RLSND173-2.5	39		
MRW2415	21	RLSND174-1.5	39		
MRW2515	21	RLSND174-2.5	39		
MRW3015	21	RLSND314-2.5	38		
NCV115	43	RLSP310	39		
NCVB100	33	RLSP350	38		
NCVD80L	34	RLSP351	38		
NCVG40L-1.0	34	RLSP430	40		
NCVG40L-2.0	34	RLSP450	38		
NCVG50	33	RLSP482	40		
NCVKIT15	33	RLSP483	40		
NCVKIT60L	34	RLSP484	40		
NCVR100	33	RLSP486	41		
NCVR40L	34	RLSP487	44		
NCVT100	33	RLSP492	43		
NCVT80L	34	RLSP503	39		
PDECUS15026	43	RLSP504	39		
PDECUS25026	43	RLSP511	30		
PDECUL15026	43	RLSP512	40		
PDECUL25026	43	RLSP513	40		
PSTCUS15026	43	RLSP516	21		
PSTCUS25026	43	RLSP518	21		
PSTCUS30026	43	RLSP527	40		
PSTCUL15026	43	RLSP531C	39		
PSTCUL25026	43	RLSP536	43		
RLBTP100-2.5	36	RLSP537	40		
RLBTP101-2.5	36	RLSPGSP120	43		
RLBTP150-2.5	36	SCI1026	26		
RLBTP151-2.5	36	SCI420	26		
RLBTP200-2.5	36	SC2315T	26		
RLBTP201-2.5	36	SC2325	26		
RLIT131-1.0	25	SC2410	26		
RLITD2910	25	SC2515T	26		
RLITCTGNM1.0	25, 42	SC2520	26		
RLITCTGNL1.0	25, 42	SC2615	26		
RLITN1515	25	SC3015	26		
RLMSP050-2.5	36	SCDECUL2225	30		
RLMSP051-2.5	36	SNB1710	40		
RLMSP100-2.5	36	SPDC1810	27		
RLMSP101-2.5	36	STCUS15026	43		
RLSND107-1.5	38	STCUS25026	43		
RLSND107-2.5	38	STCUS30026	43		
RLSND108-2.5	38	STCUL15026	43		
RLSND109-1.5	38	STCUL25026	43		
RLSND109-2.5	38	STDTL	25		
RLSND110-1.0	38	STDTM	25		
RLSND110-1.5	38	STDT5	25		

TERMS AND CONDITIONS

Ordering

We make every effort to process and ship orders as quickly as possible. Orders received by 3pm EST are typically processed the same day but please be aware that during peak ordering seasons processing times may be longer.

We do our best to have product available. Due to order fluctuations, back-orders can occur despite our best planning efforts. Please provide a forecast for your product needs to assist us with proper manufacturing planning and inventory including strategies if possible. Forecast information increases or decreases in product needs can be forwarded to our team at sales@rhythmink.com at any time. Thank you for your business and partnership as we work together to improve patient care.

To place an order,

- Call 866.633.3754 [803.252.1222]
- Fax 803.227.1015
- Email sales@rhythmink.com

Rhythmink products are currently available for sale in the United States, Australia, Brazil, Canada, Israel, Mexico and European Union [EU]. Contact our Sales Team, sales@rhythmink.com, for specific country inquiries about product availability.

Requesting Samples

Contact our sales team to request samples. About two weeks after you've received your samples, we may send you a survey to fill out that will help us to continually improve our products and ordering processes. Contact our Sales Team for details. *Not all products are part of the sample program.*

Shipping

We ship all orders via FedEx or UPS. Please let us know if you have a carrier preference at the time of your order.

Both FedEx and UPS require the following information to ensure that orders are received in a timely manner:

- Phone Number
- Address classification [residential or commercial]
- Note: FedEx and UPS do not ship to PO Boxes

Returns

No returns will be accepted without pre-authorization prior to any return and may be subject to a stocking fee. For questions regarding a return please contact us at sales@rhythmink.com.

Offer and Acceptance

Any order placed by Customer is expressly limited to and governed by these Terms and Conditions ("Agreement"). Customer's submission of an order to Rhythmink or Rhythmink's shipment of products or performance of services shall constitute Customer's agreement to and acceptance of these Terms and Conditions.

Purchase Orders

Any terms stated in, or otherwise provided in connection with, a Customer purchase order ("Purchase Order") that are in addition to or conflict with this Agreement shall have no force or effect.

Instructions for Use

Rhythmink will provide to the Customer for Customer's internal use one (1) copy of the current Instructions for Use ("IFU") for each Product. Customer shall not reproduce any IFU or Service Manual.

Delivery

Unless otherwise stated on the front of the invoice, Buyer shall be responsible for all freight, insurance, and other costs related to transport of the Products from Rhythmink's warehouse, and Buyer shall bear all risk of loss or damage from that point. Buyer shall promptly notify Rhythmink, in no event later than five (5) business days after delivery, of any claimed shortages or rejection as to any delivery. Such notice shall be in writing and shall be reasonably detailed, stating the grounds for any such rejection. Failure to give such notice within such time shall be deemed an acceptance in full of any delivery. Rhythmink shall not be liable for any shipment delays beyond the reasonable control of Rhythmink which affect Rhythmink or any of Rhythmink's suppliers, including, but not limited to, delays caused by natural disasters, acts of war, acts or omissions of Buyer, fire, strike, riot, governmental interference, unavailability or shortage of materials, labor, fuel, or power through normal commercial channels at customary and reasonable rates, or transport failures.

TERMS AND CONDITIONS

Price and Payment

RhythmLink, at its discretion, may require reasonable advance assurances of payment. All unpaid invoices shall bear interest at an amount equal to 2% of the outstanding balance per month (or the maximum rate of interest allowed to be contracted for by law, whichever is less) commencing upon the date payment is due. Buyer's failure to make timely payment may result in such action as commencement of proceedings for collection, revocation of credit, stoppage of shipment, delay or cessation of future deliveries, repossession of unpaid delivered goods, and termination of any one or more sales agreements. RhythmLink retains (and Buyer grants to RhythmLink by submitting a purchase order) a security interest in the Products to secure payment in full and compliance with all sales agreements, and Buyer agrees to execute any additional documents necessary to perfect such security interest and authorizes RhythmLink to execute and file financing statements or other necessary documents to perfect and continue such security interest. In the event that a sales invoice shall be placed by RhythmLink in the hands of an attorney for the purpose of collection or for the purpose of enforcing RhythmLink's security interest in the Products, Buyer agrees to pay any and all costs associated with such placement, including without limitation, attorney's fees and costs incurred prior to, during, or subsequent to trial. The amount due RhythmLink is contracted in U.S. Dollars and payment in U.S. Dollars is of the essence.

Taxes

Prices and fees stated for Products do not include any taxes. Customer shall pay when due any applicable sales or other taxes and related interest and penalties arising from Customer's acquisition or possession of the Products.

Restriction on Use of Products

Customer shall use each Product only: (i) for Customer's internal use; (ii) in the manner described in the Product's IFU; (iii) in accordance with applicable laws and regulations, and (iv) shall not export or re-export technical data or Products in violation of applicable export regulations.

Warranty

RhythmLink warrants that the products shall be new and of good quality and free from defects in material and workmanship for a period of 12 months from the date of delivery. In the event any of the products are found to be defective during such., RhythmLink will, at its option, repair or replace the defective product or refund the purchase price attributable. If a product return is required a Return Material Authorization (RMA) will be issued and the product must be returned within 90 days of the RMA being issued or the authorization to return could be revoked.

EXCEPT AS PROVIDED ABOVE, RHYTHMLINK MAKES NO REPRESENTATION OR WARRANTIES OF ANY KIND, EXPRESS OR IMPLIED, WITH RESPECT TO THE CONDITION OR PERFORMANCE OF THE PRODUCTS AND IT HEREBY DISCLAIMS ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. UNDER NO CIRCUMSTANCES SHALL RHYTHMLINK BE LIABLE FOR ANY INDIRECT, INCIDENTAL, CONSEQUENTIAL, SPECIAL OR PUNITIVE LOSSES OR DAMAGES, INCLUDING THOSE RESULTING FROM THE PRODUCTS, THEIR DELIVERY, OR THEIR USE.

Limitation of Liability

RhythmLink shall not be liable under any circumstances for any special, consequential, incidental, punitive or exemplary damages arising out of or in any way connected with the products, including, but not limited to, damages for lost profits, loss of use, or for any damages or sums paid by buyer to third parties, even if RhythmLink has been advised of the possibility of such damages.

CAUTION!! Federal Law restricts this device to sale by or on the order of a practitioner licensed by the Law of the State in which he/she practices.

Quality

Have confidence in the products that you purchase from RhythmLink. We stand behind all of our products and have earned some of the most important certifications from around the globe.

- ISO 13485:2016
- EN ISO 13485:2016
- MDD 93/42/EEC
- MDSAP (ISO 13485:2016)
- ISO Class 8 Cleanroom

Contact Quality or Sales for specific marketing authorization.

For more information visit RhythmLink.com/quality.

Trademarks

RhythmLink is a registered trademark of RhythmLink International, LLC.

Quick Connect System, Sticky Pad, Limelight, LeadLock, Slim Cup, Deep Cup and Webb are trademarks of RhythmLink International, LLC.

The appearance of the MR Quick Connect extension cable is a trademark of RhythmLink.

Ten20 and Nuprep are registered trademarks of Weaver and Co.

Dragonfly and Cobra are registered trademarks of Neurovision Medical Products.

SightSaver is a trademark of AnsheI Technology Inc.

Patents

PressOn Electrode, US Patent Nos. 8112131, 8185183, 8271065, 8428681, 8538502, Chinese Patent No. CN101720202, Japanese Patent No. JP5226776 and European Patent No. EP2155056.

Monopolar Stimulating Probe with Removable Handle, US Patent No. 8137285

Guardian Needle Electrode, US Patent Nos. 10,932,683 and 11,690,549

Enhancing Patient Care Worldwide by
Connecting Patients to Machines.



PHONE
866.633.3754
[803.252.1222]

FIND US ONLINE AT
RhythmLink.com
sales@rhythmLink.com

